

Product datasheet for **SC310408**

RBM34 (NM_015014) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RBM34 (NM_015014) Human Untagged Clone
Tag:	Tag Free
Symbol:	RBM34
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC310408 representing NM_015014. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCTTGAAGGGATGAGCAAACGGAAGAGAAAGAGAAGTGTCCAGGAGGGAGAGAATCCTGACGAC
GGCGTTCGCGGGAGTCCGCCGGAAGACTACAGGCTTGGACAGGTCCGCGAGTAGCTTATTTCCGCGCGAA
CACCATTCCAGAGGTGGCACCCTGGCTGGCGTCCCTCTTCAGTTCTCTGGAGCCCCAGATTCAACCC
GTGTACGTGCCTGTGCCTAAACAACCATCAAAAAACGAAACGGAATGAGGAGGAAGAAAGTACATCC
CAGATTTAAAGACCCTTTCCGAAGAACCTGCCAAAAAGTGAAGCGAAGAAGAAACACACTAACGCA
GAAAAAAGTTGGCAGACAGGAAAGCGCTCTAGCGAGTCTGATTTAGAAGAAGAAATTCACAGAAA
CAAGGGCAGAAAAGGAAAAATTCTCAACCTGGTGTAAAGTAGCAGATAGAAAAATACTTGATGACACA
GAAGACACAGTTGTCAAGTCAAAGAAAGAAAATTCAAATCAACCAAGAAAGAGAGATTAAGAATGAG
AGAAGTGTGTTTGGGAATTTGCCTGTACATGTAATAAGAAGAAGCTGAAGTCGTTTTTAAAGAG
TATGGACAAATAGAATCTGTACGATTTCTGTTCTCTGATTCCAGCAGAGGGAACGCTATCCAAAAAGTTG
GCAGCAATAAACGTAATTCATCCTGATCAGAAAAATATTAATGCCTATGTTGTGTTTAAAGGAGGAG
AGTGTGCCACGCAAGCATTGAAAAGAAATGGGGCCAGATTGCAGATGGATTTTCGATTAGAGTTGAT
CTCGCATCTGAGACCTCATCTAGAGACAAGAGATCGGTTTTTGTGGGAATCTCCCTATAAAGTTGAA
GAATCTGCCATTGAGAAGCACTTCTGGACTGTGGAAGTATCATGGCCGTGAGGATTGTGAGAGACAAA
ATGACAGGCATCGCAAAGGGTTTTGGCTATGTGCTCTTTGAGAATACAGATTCTGTTTCATCTTGCTCTG
AAATTAATAATTCTGAACTCATGGGGAGAAAACCTCAGAGTCATGCGTTCTGTTAATAAAGAAAAATTT
AAACAACAAAATCAAATCCACGATTGAAGAATGTCAAGAAACCTAAGCAGGGACTTAATTTTACTTCC
AAAAGTGCAGAAGGACATCTAAAGCTTATTTATTGGAGAAAAAGCTGTTCTCCTTAAACGAAGAAG
AAAGGACAGAAGAAAAGTGGACGCCCTAAGAAACAGAGAAAACAGAATAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online >](#)

Plasmid Map:	□
ACCN:	NM_015014
Insert Size:	1293 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_015014.2
RefSeq Size:	1870 bp
RefSeq ORF:	1293 bp
Locus ID:	23029
UniProt ID:	P42696
Cytogenetics:	1q42.3
MW:	48.6 kDa
Gene Summary:	<p>This gene encodes a member of the RNA-binding motif family of RNA recognition motif proteins. The encoded protein contains an RNA-binding domain made up of two RNA recognition motif subdomains referred to as RNA recognition motif-1 and RNA recognition motif-2. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2016]</p> <p>Transcript Variant: This variant (1) encodes the longest isoform (1).</p>