

Product datasheet for **SC310416**

FDPS (NM_002004) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: FDPS (NM_002004) Human Untagged Clone
Tag: Tag Free
Symbol: FDPS
Synonyms: FPPS; FPS; POROK9
Mammalian Cell Selection: None
Vector: pCMV6-XL5
E. coli Selection: Ampicillin (100 ug/mL)

Fully Sequenced ORF: >OriGene sequence for NM_002004 edited
 ATGCCCTGTCCCCTGGTTGAGATCTGTGGGGTCTTCCTGCTGCCAGCCCCCTACTGG
 GCACCCCGGAGAGGTGGCTGGTTCCCTACGGCGGCCCTCCCTGGTGCACGGGTACCCA
 GTCCTGGCCTGGCACAGTGCCCGCTGCTGGTGCCAAGCGTGACAGAGGAACCTCGAGCC
 CTTTGCTCCTCCCTCAGAATGAACGGAGACCAGAATTCAGATGTTTATGCCCAAGAAAAG
 CAGGATTCGTTCACTTCTCCAGATCGTTAGGGTGTGACTGAGGATGAGATGGGG
 CACCCAGAGATAGGAGATGCTATTGCCCGGCTCAAGGAGGTCTGGAGTACAATGCCATT
 GGAGGCAAGTATAACCGGGTTTGACGGTGGTAGTAGCATTCCGGGAGCTGGTGGAGCCA
 AGGAAACAGGATGCTGATAGTCTCCAGCGGCCCTGGACTGTGGGCTGGTGTGGAACTG
 CTGCAAGCTTTCTCCTGGTGGCAGATGACATCATGGATTCATCCCTTACCCGCCGGGA
 CAGATCTGCTGGTATCAGAAGCCGGCGTGGGTTGGATGCCATCAATGATGCTAACCTC
 CTGGAAGCATGTATCTACCGCTGCTGAAGCTCTATTGCCGGGAGCAGCCCTATTACCTG
 AACCTGATCGAGCTTCTCCTGCAGAGTTCCTATCAGACTGAGATTGGGCAGACCCCTGGAC
 CTCCTCACAGCCCCCAGGGCAATGTGGATCTTGTGAGATTCACTGAAAAGAGGTACAAA
 TCTATTGTCAAGTACAAGACAGCTTTCTACTCCTTCTACCTTCTATAGCTGCAGCCATG
 TACATGGCAGGAATTGATGGCGAGAAGGAGCACGCCAATGCCAAGAAGATCCTGCTGGAG
 ATGGGGGAGTTCTTTTCAAGTTCAGGATGATTACCTTGACCTCTTTGGGGACCCAGTGTG
 ACCGGCAAAATTGGCACTGACATCCAGGACAACAAATGCAGCTGGCTGGTGGTTCAGTGT
 CTGCAACGGGCCACTCCAGAACAGTACCAGATCCTGAAGGAAAATTACGGGCAGAAGGAG
 GCTGAGAAAAGTGCCCGGTGAAGGCGCTATATGAGGAGCTGGATCTGCCAGCAGTGTTCT
 TTGCAATATGAGGAAGACAGTTACAGCCACATTATGGCTCTCATTGAACAGTACGCAGCA
 CCCCTGCCCCAGCCGTCTTCTGGGGCTTGCGCGCAAAATCTACAAGCGGAGAAAGTGA

Restriction Sites: Please inquire
ACCN: NM_002004
Insert Size: 1300 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_002004.2 , NP_001995.1
RefSeq Size:	1345 bp
RefSeq ORF:	1260 bp
Locus ID:	2224
UniProt ID:	P14324
Cytogenetics:	1q22
Domains:	polyprenyl_synt
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Terpenoid backbone biosynthesis
Gene Summary:	<p>This gene encodes an enzyme that catalyzes the production of geranyl pyrophosphate and farnesyl pyrophosphate from isopentenyl pyrophosphate and dimethylallyl pyrophosphate. The resulting product, farnesyl pyrophosphate, is a key intermediate in cholesterol and sterol biosynthesis, a substrate for protein farnesylation and geranylgeranylation, and a ligand or agonist for certain hormone receptors and growth receptors. Drugs that inhibit this enzyme prevent the post-translational modifications of small GTPases and have been used to treat diseases related to bone resorption. Multiple pseudogenes have been found on chromosomes 1, 7, 14, 15, 21 and X. Multiple transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Oct 2008]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longer isoform (a). Both variants 1 and 2 encode the same isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.</p>