

Product datasheet for **SC310417**

RBFOX1 (NM_145891) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RBFOX1 (NM_145891) Human Untagged Clone
Tag:	Tag Free
Symbol:	RBFOX1
Synonyms:	2BP1; A2BP1; FOX-1; FOX1; HRNBP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC310417 representing NM_145891. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTGGCGTCTCAAGGAGTTCTCCTGCATCCTTATGGCGTGCCTATGATTGTACCGGCAGCTCCTTAC
CTTCTGGACTGATTCAGGGTAATCAGGAAGCAGCCGCTGCCCTGACACAATGGCTCAGCCTTACGCT
TCGGCCAGTTTCTCCCCGAGAACGGTATCCCCGCGAATACACGGCCCTCATCCCCACCCCGCG
CCAGAGTACACAGGCCAGACCACGGTCCCGAGCACACATTAACCTGTACCCTCCGCCCAGACGCAC
TCCGAGCAGAGCCCGCGGACACGAGCGCTCAGACCGTCTTGGCACCGCCACACAGACAGATGACGCA
GCACCGACGGATGCCAGCCCCAGACACAACCTTCTGAAAACACGAAAACAAGTCTCAGCCCAAGCGG
CTGCATGTCTCAATATCCCCTTCCAGTTCCGGGATCCGGACCTCAGACAATGTTTGGTCAATTTGGT
AAAATCTTAGATGTTGAAATTTTTTAAATGAGCGAGGCTCAAAGGGATTTGGTTTCGTAACCTTCGAA
AATAGTGCCGATGCGGACAGGGCGAGGAGAAATTACACGGCACCGTGGTAGAGGGCCGTAATTCGAG
GTAATAATGCCACAGCACGTGTAATGACAAAATAAAAAGACCGTCAACCTTATACAAAATGGCTGGAAA
TTGAATCCAGTTGTGGGTGCAGTCTACAGTCCCGAATTCTATGCAGGCACGGTCTGTGTGCCAGGCC
AACCAGGAGGGATCTCCATGTACAGTGCCCCAGTTCCTTGTATATACTTCTGCAATGCCAGGCTTC
CCGATCCAGCAGCCACCGCCGCGGCCCTACCGAGGGGCGCACCTGCGAGGCCGCGGTGCGACCGTG
TACAACACCTTCAGGGCCGCGGCCGCCCGCCCGATCCCGGCTACGGCGGAGTAGTGTATCAAGAG
CCTGTGTATGGCAATAAATTGCTGCAGGGTGGTTATGCTGCATACCGCTACGCCACGCTACCCCTGCC
ACTGCCGTGCCTACAGTGACAGTTACGGACGAGTTTATGCTGCCGACCCCTACCACCACGCACTTGCT
CCAGCCCCACCTACGGCGTTGGTGCCATGAATGCTTTTGCACCTTTGACTGATGCCAAGACTAGGAGC
CATGCTGATGATGGGTCTCGTTCTTTCTTATTGCAGGCTAGTATATACCGAGGGGATACAACCGT
TTTGCTCCATACTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
Plasmid Map:	□
ACCN:	NM_145891
Insert Size:	1257 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_145891.2
RefSeq Size:	4102 bp
RefSeq ORF:	1257 bp
Locus ID:	54715
UniProt ID:	Q9NWB1
Cytogenetics:	16p13.3
Domains:	RRM
MW:	44.8 kDa

Gene Summary:

The Fox-1 family of RNA-binding proteins is evolutionarily conserved, and regulates tissue-specific alternative splicing in metazoa. Fox-1 recognizes a (U)GCAUG stretch in regulated exons or in flanking introns. The protein binds to the C-terminus of ataxin-2 and may contribute to the restricted pathology of spinocerebellar ataxia type 2 (SCA2). Ataxin-2 is the product of the SCA2 gene which causes familial neurodegenerative diseases. Fox-1 and ataxin-2 are both localized in the trans-Golgi network. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2011]

Transcript Variant: This variant (1), also known as gamma, encodes isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data because no single transcript was available for the full length of the gene. The extent of this transcript is supported by transcript alignments.