

Product datasheet for SC310568

IFI35 (NM_005533) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	IFI35 (NM_005533) Human Untagged Clone
Tag:	Tag Free
Symbol:	IFI35
Synonyms:	IFP35
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC310568 representing NM_005533. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCAGCCCCACTGGATGCCGCCCTCCACGCCCTTCAGGAGGAGCAGGCCAGACTCAAGATGAGGCTG
TGGGACCTGCAGCAGCTGAGAAAGGAGCTCGGGGACTCCCCAAAGACAAGGTCCCATTTTCAGTGCCC
AAGATCCCCCTGGTATTCCGAGGACACACCCAGCAGGACCCGGAAGTGCCTAAGTCTTTAGTTTCCAAT
TTGCGGATCCACTGCCCTCTGCTTGGGGCTCTGCTCTGATCACCTTTGATGACCCCAAAGTGGCTGAG
CAGGTGCTGCAACAAAAGGAGCACACGATCAACATGGAGGAGTGCCGGCTGCGGGTGCAGGTCCAGCCC
TTGGAGCTGCCCATGGTACCACCATCCAGGTGATGATGTCCAGCCAGTTGAGTGCCCGGAGGGTGTG
GTCACTGGATTTCTGCCAGCCTCAGGCTGAGTGAGGAGGAGCTGCTGGACAAGCTAGAGATCTTCTTT
GGCAAGACTAGGAACGGAGGTGGCGATGTGGACGTTTCGGGAGCTACTGCCAGGGAGTGTCATGCTGGGG
TTTGCTAGGGATGGAGTGGCTCAGCGTCTGTGCCAAATCGGCCAGTTCACAGTGCCACTGGTGGGCAG
CAAGTCCCCTCTGAGAGTCTCTCCGTATGTGAACGGGGAGATCCAGAAGGCTGAGATCAGGTGCGAGCCA
GTTCCCGCTCGGTACTGGTGTCAACATTCCTGATATCTTGGATGGCCCGGAGCTGCATGACGTCTCTG
GAGATCCACTTCCAGAAGCCACCCGCGGGGGCGGGGAGGTAGAGGCCCTGACAGTCGTACCCCAAGGA
CAGCAGGGCCTAGCAGTCTTACCTCTGAGTCAGGCTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_005533

Insert Size: 867 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_005533.4](#)

RefSeq Size: 1260 bp

RefSeq ORF: 867 bp

Locus ID: 3430
UniProt ID: [P80217](#)
Cytogenetics: 17q21.31
MW: 31.8 kDa
Gene Summary: Not yet known.[UniProtKB/Swiss-Prot Function]