

Product datasheet for **SC310596**

CCDC90B (NM_021825) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	CCDC90B (NM_021825) Human Untagged Clone
Tag:	Tag Free
Symbol:	CCDC90B
Synonyms:	MDS011; MDS025
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC310596 representing NM_021825. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAATAGTCGCCAGGCTTGGCGGCTCTTTCTCTCCCAAGGCAGAGGAGATCGTTGGGTTTCAAGGCC
CGCGGGCATTCTCGCCGGCCCTGCGGAGAGATTCTTCACTACCACAACCAAGGAGGATATGATAGG
CGGCCAGTGGATATAACTCCTTTAGAACAAGGAAATTAACCTTTGATACCCATGCATTGGTTCAAGGAC
TTGGAACTCATGGATTTGACAAAACACAAGCAGAAACAATTGTATCAGCGTTAACTGCTTTATCAAAT
GTCAGCCTGGATACTATCTATAAAGAGATGGTCACTCAAGCTCAACAGGAAATAACAGTACAACAGCTA
ATGGCTCATTGGATGCTATCAGGAAAGACATGGTCATCCTAGAGAAAAGTGAATTTGCAATCTGAGA
GCAGAGAATGAGAAAATGAAAATTGAATTAGACCAAGTTAAGCAACAATAATGCATGAAACCAGTCGA
ATCAGAGCAGATAATAAACTGGATATCAACTTAGAAAAGGAGCAGAGTAACAGATATGTTTACAGATCAA
GAAAAGCAACTTATGAAAACAACACAGAAATTTACAAAAAGGATACTCAAACAAAAGTATTATTTCA
GAGACCAGTAATAAAATGACGCTGAAATTGCTTCCTTAAAAACTGATGGAATCTAACAACTTGAG
ACAATTCGTTATCTTGACGCTTCGGTGTACTTGGCTGGCAATAGCATTGGGATTTATAGATTCTGG
AAGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_021825

Insert Size: 765 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_021825.4](#)

RefSeq Size: 4018 bp

RefSeq ORF: 765 bp
Locus ID: 60492
UniProt ID: Q9GZT6
Cytogenetics: 11q14.1
Protein Families: Transmembrane
MW: 29.5 kDa