

Product datasheet for SC310702

RNF151 (NM_174903) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RNF151 (NM_174903) Human Untagged Clone
Tag:	Tag Free
Symbol:	RNF151
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC310702 representing NM_174903. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
 ATGGGTGGCGGTATGATCTCAACCTCTTCGCCAGCCCTCCTGACAGCAACTTCGTGTGCTCCGTCTGC
 CATGGGTTTCTCAAGAGGCCAGCAAGTTGCCATGCAGCCACATCTTCTGCAAAAAGTGCATCCTCCGG
 TGGCTAGCCAGACAAAAGACCTGTCCGTGCTGTAGGAAAGAGGTGAAAAGGAAAAGTTGTCCACATG
 AATAAACTCCGAAAACCATTTGGCCGCTGGAAGTCAAGTGAAGAACGCCGACGCTGGCTGCATAGTG
 ACATGCCCCCTGGCCCATCGAAGGGGCACCAGGACTCATGCCCTTTGAGCTAACGGCCTGCCCAAC
 GAGGGCTGCACCTCGCAGGTGCCGCTGGGACCCTGGCAGAGCACCAGGAGCATTGCCAGCAAGGGTCC
 CAGCAGCGCTGCCCTTGGCTGCGGGGCCACCCTGGACCCGGCCGAGCGTGCTCGCCACAACCTGTAC
 CGGGAGCTGCACAACGCCTGGAGCGTGCGCCAGGAGCGCCGTCGGCCCTGCTGTGTCCCTCCTGCGG
 CGTGTGCGCTGGCTGGACCAAGCCACCAGTGTCTTCGTAGAGAGCTGGCGGAGCTCAGCAACTCCTG
 GAGGAAGACACCGCTCTGCTGGAGGGTGGCCACAGGAGGAGCCGAGGCTGCCCCAGAAGGCAACGTT
 GGGGCTGAGGTGGTGGGGAGCCAGGGCCAACATACCTTGTAATAG
 ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	Sgfl-MluI
ACCN:	NM_174903
Insert Size:	738 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_174903.4</u>
RefSeq Size:	796 bp
RefSeq ORF:	738 bp
Locus ID:	146310
UniProt ID:	<u>Q2KHN1</u>
Cytogenetics:	16p13.3
Protein Families:	Druggable Genome
MW:	27.4 kDa
Gene Summary:	May be involved in acrosome formation of spermatids.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longer isoform (1).