

Product datasheet for **SC310744**

ZNF705D (NM_001039615) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF705D (NM_001039615) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF705D
Synonyms:	ZNF705C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC310744 representing NM_001039615. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
 ATGCATTCAGTAGAGAAAGTGACTTTTGAAGATGTAGCTATTGACTTCACCCAGGAAGAGTGGGACATG
 ATGGACACATCCAAGAGAAAGCTGTACAGAGATGTGATGCTGGAAAAATCAGTCACCTGGTGTCCCTC
 GGGTACCAGATAAGCAAATCCTATATAATTTTGCAGCTGGAGCAAGGAAAAGAGCTGTGGCGGGAAGGA
 AGAGTATTTCTTCAAGACCAGATCCAGACAGGAAAGTGCCCTTAAGAAAAACACATGATATCCATG
 CATCCTATCATCAGAAAAGACGCATCCACAGTATGACAATGGAGAACTCTCTATTCTGGAGGATCCT
 TTTGAATATAATGATTCGGGAGAAGATTGCACTCACAGTTCACAATAACTCAGTGTGTTGTTAACTCAC
 AGTGGAAGAAACCCTGTGTGAGCAAAACAGTGTGGAAAAATCCCTTCGTAATCTTTTGTCCCTAAACCA
 CGTAAACAAATTCATACTAAAGGTAAATCATATCAATGTAATCTATGTGAAAAGGCCTATACTAATTGC
 TTTTACCTTAGACGGCACAAGATGACTCACACTGGAGAGAGGCCATATGCATGTCATCTATGTGAAAAA
 GCCTTCACTCAGTGTTCTCACCTTAGAAGACATGAGAAAACACACGGGAGAGAGACCATATAAGTGT
 CATCAATGTGGAAAGCCTTTATTCAATCCTTTAACCTTCGAAGACATGAGAGAACTCACCTTGGACAA
 AAGTGTATGAATGTGATAAAAGTGGGAAAGCCTTTAGTCAAAGCTCTGGCTTTAGAGGAAACAAAATA
 ATTCACATTGGAGAGAAACCACATGCTTGTCTTCTATGTGGGAAGGCCTTCAGTCTGTCTCCGACCTT
 AGATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001039615
Insert Size:	903 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001039615.3
RefSeq Size:	3683 bp
RefSeq ORF:	903 bp
Locus ID:	728957
UniProt ID:	P0CH99
Cytogenetics:	8p23.1
MW:	34.7 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]