

Product datasheet for **SC310767**

MAP2 (NM_001039538) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MAP2 (NM_001039538) Human Untagged Clone
Tag:	Tag Free
Symbol:	MAP2
Synonyms:	MAP-2; MAP2A; MAP2B; MAP2C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC310767 representing NM_001039538.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGCAGATGAACGGAAAGATGAAGCAAAGGCACCTCACTGGACCTCAGCACCGCTAACAGAGGCATCT
GCACACTCACATCCACCTGAGATTAAGGATCAAGGCGGAGCAGGGGAAGGACTTGTCCGAAGCGCCAAT
GGATTCCCATACAGGGAGGATGAAGAGGGTGCCTTTGGAGAGCATGGGTACAGGGCACCTATTCAAAT
ACCAAAGAGAATGGGATCAACGGAGAGCTGACCTCAGCTGACAGAGAAACAGCAGAGGAGGTGTCTGCA
AGGATAGTTCAAGTAGTCACTGTGAGGCTGTAGCAGTCTGAAAGGTGAACAAGAGAAAGAAGCTCAA
CATAAAGACCAGACTGCAGCTCTGCCTTAGCAGCTGAAGAAACAGCTAATCTGCCTCCTTCTCCACCC
CCATCACCTGCCTCAGAACAGACTGTACAGTGGAGGAAGCAGCAGGTGGGGAATCAGCTCTGGCTCCC
AGTGTATTTAAACAGGCAAAGGACAAAGTCTCTAATTCTACCTTGTCAAAGATTCCTGCTTACAGGGT
AGCACAAAGTCCCAAGATACAGCTCAGCTGCCCTAGCAGACTAAAAGGGCTACATTTTCTGACAGT
TTATTAATACAGCCACCTCAGCAGGCTCCACAGACCGTTTGCCATACTCAAATCAGGGAAACAAGGAC
GGAGTAACCAAGAGCCCAGAAAAGCGCTTCTCTCCCAAGACCTTCTCCATTCTCCCTCCTCGGCCGA
GGTGTGTCAAGGACAGAGATGAGAATTCCTTCTCTCAACAGTTCTATCTCTTCTCAGCACGGCGG
ACCACCAGGTCAAGCCAATTCGAGAGCAGGGAAGAGTGGTACCTCAACACCCACTACCCTGGGTCT
ACTGCCATCACTCCTGGCACCCCAAGTTATTTCTTACGCACACCAGGCACTCCTGGAACCCCTAGC
TATCCCAGGACCCTCACACACCAGGAACCCCAAGTCTGCCATCTTGGTGCCGAGTGAGAAGAAGGTC
GCCATCATACGTACTCTCCAAAATCTCCTGCGACTCCCAAGCAGCTTCGGCTTATTAACCAACCACTG
CCAGACCTGAAGAATGTCAAATCCAAAATCGGATCAACAGACAACATCAAATACCAGCCTAAAGGGGGG
CAGGTTAGGATTTTAAACAAGAAGATCGATTTTAGCAAAGTTCAGTCCAGATGTGGTTCCAAGGATAAC
ATCAAACATTTCGGCTGGGGCGGAAATGTACAAATTGTTACCAAGAAAATAGACCTAAGCCATGTGACA
TCCAAATGTGGCTCTCTGAAGAACATCCGCCACAGGCCAGGTGGCGGACGTGTGAAAATTGAGAGTGTA
AACTAGATTTCAAAGAAAAGGCCAAGCTAAAGTTGGTTCTCTTGATAATGCTCATCATGTACCTGGA
GGTGGTAAATGTCAAGATTGACAGCCAAAAGTTGAACTTCAGAGAGCATGCTAAAGCCCGTGTGGACCAT
GGGGCTGAGATCATTACACAGTCCCAGGCAGATCCAGCGTGGCATCACCCGACGACTCAGCAATGTC
TCCTCGTCTGGAAGCATCAACTGCTCGAATCTCCTCAGCTTGCCACTTTGGCTGAGGATGCTACTGCT
GCACTCGCTAAGCAGGGCTTTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001039538

Insert Size: 1680 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001039538.1](#)

RefSeq Size: 5844 bp

RefSeq ORF: 1680 bp

Locus ID: 4133

UniProt ID: [P11137](#)

Cytogenetics: 2q34

Protein Families:	Adult stem cells, Druggable Genome, ES Cell Differentiation/IPS
MW:	59 kDa
Gene Summary:	<p>This gene encodes a protein that belongs to the microtubule-associated protein family. The proteins of this family are thought to be involved in microtubule assembly, which is an essential step in neurogenesis. The products of similar genes in rat and mouse are neuron-specific cytoskeletal proteins that are enriched in dendrites, implicating a role in determining and stabilizing dendritic shape during neuron development. A number of alternatively spliced variants encoding distinct isoforms have been described. [provided by RefSeq, Jan 2010]</p> <p>Transcript Variant: This variant (5) differs in the 5' UTR and includes an alternate in-frame exon, compared to variant 2, resulting in a longer protein (isoform 5), compared to isoform 2.</p>