

Product datasheet for **SC310795**

TMEM232 (NM_001039763) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM232 (NM_001039763) Human Untagged Clone
Tag:	Tag Free
Symbol:	TMEM232
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC310795 representing NM_001039763.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAATATGCCTGTTAACAAATCACCTATGATTAATACATGTGGAGGCATATCTCCCTTATCATGAA
GAGCTCTGGAAATTAATTTTCAACATTTAAGTGGAGAAAGGGTCATAAATCAAGGCCAACATTTTCA
ATCACAAAAGAATTCATCTTGAGATTCATCAAAACACAAAATTCAAAGAGAAGGAAGAATTGTTGGAA
CTAGCTAGAAAAATCATTCTCAGATGTAAGAGAAAATTGGGTCTCAAAACCTGGGCTCTGGAAGGCAT
GTGCATCTTCTGCTGCATGGACTGAAGTAATATATCTGGCTCAATGCAAAGGGGAAATCCAAGATGAA
TCTTTAAATATGCTTTATGCATCTCTGGACCATGCTTCTTTGATTATGATCATCTGCCTGCCTTATTT
TTTGTGCGGAATCGGTTTTATACAGACTGTGTTGTGATGCATCCCTCAAAACATATTTATATTAGTT
GAAATAAGCTAGCAAAGATTGGCTATTTGGTCTTCTACGACTTTTTATTTTTCTGCATGGTCAT
CTAGAAAGTTTTAAACAACATTTACTTAGGCTTCAACCATATTTATACGCACTTTCTTTCTGAGCA
TCATATCACAGTATCCAAACATCTTTCAAATGTGAATTCATCCTGAAAGCCTCGGAAATTATAGGT
AAAAGAGAACTCCGTTCTGAATCCATTTTAGACCTGTGGAAGATAAGAAAAGATATGAGAACACAGAT
TCTGATATGGGAGGATATGAAATTAACCACCTGCTCTGGCACTGTGTTGCTGCTTGGTCTTGTGTTGAG
AATAACAGTCTCAGTTGAATAACGTGCTTGAACATCTCGTCTCCATAAGACACAGCTTCAAAAGAAA
TGCTGGTTGGATTAGTACTGGCTTTACTGGTCTTGGGAGGCTGCCAAATTAACATGGCCTGCTTG
AAAGCTTTAATGGACGTAGTGAGAGATTTGTTTCAAGCATTATGTCTGTTCAAAATCAGGAAGAAAGT
TGCAAGGTAGATGATTTTCTGGGCTTGAATGTAGTCTACATATATACAGTAATCTTGCAGAAATC
TGCTTGTATGCAGCCACTTCTGATTTGCGAAAACTGCTTAAATTGGTTTCTGTCACTGTAAAGTTCA
CAAAAAAATATTTTACTTGGACAAATCAGTACCTCCAGAATTAAGGAAACAAGTATTTTAAGTCTT
TTGGAATATTTCTCTTCAAAAATGTGAGAGAACTGTGATCAAGTAGTCTGGACTGGTTACTATGGCTTA
GTGTATAACCTGGTGAAATTTTATGGGAACCTCAAGGAGACGAAGAACAGGATGGACTTAGAAACATG
ATATGGCAACATTGCAGAAAACAAAGGATTATGAGGAAGATGTACGAATCCAAATGCAATCAATATA
GCTCAGGCTGAGTTAAATGACCCAATGATCCTTTCACTAGATATAGTACAAATATTTTATCAAAATGTA
GGAGAAGAAGTTTCTCCAAATATATTGGGTGGAGAATTGCCAACACTCTTCCAACTATTCTTTCCC
CCCATTGAGGCCATTTTCTCTTTGAAGAAACCATCAATAAAAAAGGATCAAAACAAATATCCAAAT
AAAAAGCTGGAGTCTGTGAAAAAGCAAGTTCTACATTTCACTGTAAGGGAACATCCTTCTGTATCAGAA
ATTCCCATGTTTCCATATCCAGATTTCTCACCAAGGCAGATAAAGAGTTGGCAAAATCATTGACCAT
CACTGGCAGGAAGAGCTAAAGATCCGAGAAAAAGAGATGCAATATGCAAGGCCCAAGAAGCTTAAAGAT
AAAAAGTTAGCAGAGAAAAATCACTTTCAAGAAGTTATGAAGAAAAGAGAAGAGAACTACATAAGCAA
ACCAAGCCCTATGAACCTCCTTAGGAAGGAAGTAATTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001039763

Insert Size: 1974 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001039763.3</u>
RefSeq Size:	3193 bp
RefSeq ORF:	1974 bp
Locus ID:	642987
UniProt ID:	<u>C9JQI7</u>
Cytogenetics:	5q22.1
Protein Families:	Transmembrane
MW:	76.5 kDa