

## Product datasheet for SC310811

### ENTR1 (NM\_001039707) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ENTR1 (NM_001039707) Human Untagged Clone
Tag:	Tag Free
Symbol:	ENTR1
Synonyms:	NY-CO-3; SDCCAG3; SDDAG3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC310811 representing NM_001039707. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCGGGCTACCGAGCGCCGCCCGGGCGCCACCCCGCTGTCCCGAGCCGGAGCCTCGCCATTCCCGAC
GCTCCAGCGTTCTATGAGCGCCGGTCTTGTCTCCCCAGCTAAATTGTGAGCGCCCCATGGCAGGGAC
CTGGACTCCCCCTTCTTCGGCATTTCGGCCGGCCTTTATGTGCTATGTGCCAGCCGGTGTGGCTTCC
GTGGGAGACACAGATTTGGCTATGAAAGGGGAAATGTTCTAAGCAGAGCCCGTCAGGAGCCACGGG
ACACATTTTGGAGATGACAGATTTGAAGATCTGGAAGAGGCAAATCCATTCTCTTTAGAGATTTCTG
AAGACCAAGAACCTCGGCCTCTCGAAAGAGGATCCGGCCAGCAGAATTTATGCAAAGGAAGCCTCGAGG
CATTCCCTGGGACTTGACCACAACCTCCCCACCCTCCCAAACCGCGGGTATGGCCTGGAGTATCAGCAG
CCATTTTTCGAGGATCCGACAGGGGCTGGTGACCTCCTGGATGAGGAGGAGGATGAGGACACCGGATGG
AGTGGGGCCTACCTGCCGTCCGCATCGAGCAGACTCACCCCGAGAGGGTCCCTGCCGGCAGCTCGCCC
TGCAGCACATACTTTCTTTTCTCCACCCCGTCGGAGCTGGCAGGGCCTGAGTCTCTGCCCTCGTGG
GCGTTGAGTGACACTGATTCTCGCGTGTCTCCGGCCTCTCCGGCAGGGAGTCTAGCGCAGACTTTGCG
GTTTCATGGAGAGTCTCTGGGAGACAGGCACCTGCCGACGCTGCAGATAAGTTACGACGCACTGAAAGAT
GAAAATTTAAGCTGAGAAGAAAGCTGAATGAGGTTTCAAGCTTCTCTGAAGCTCAAACAGAAATGGTG
AGGACGCTTGAGCGGAAGTTAGAAGCAAAAATGATCAAGGAAAGCGACTACCACGACTGGAGTCCG
GTGGTTACGAGGTGGAGCAGAACTGGAGCTGATGACCAAACGGGCTGTAAAGGCAGAAAACAGGTC
GTGAAACTAAAACAGAAATCAGTTTGCTCCAGGCGCAGGTCTCCAACCTCCAGCGAGAGAAATGAAGCC
CTGCGGTGCGGCCAGGGTCCAGCCTGACCGTGGTGAAGCAGAACGCCGACGCTGGCCCTGCAGAACCTC
CGGGTGGTCATGAACAGTGCACAGGCTTCCATCAAGCAACTGGTTTCCGGAGCTGAGACTGAATCTT
GTTGCCGAAATCCTTAATCTATAGACAGAATTTCTGAAGTTAAAGACGAGGAGGAAGACTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM\_001039707

Insert Size: 1308 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001039707.1](#)

RefSeq Size: 2390 bp

RefSeq ORF: 1308 bp

Locus ID: 10807

UniProt ID: [Q96C92](#)

Cytogenetics: 9q34.3

MW: 48 kDa

**Gene Summary:**

Endosome-associated protein that plays a role in membrane receptor sorting, cytokinesis and ciliogenesis (PubMed:23108400, PubMed:25278552, PubMed:27767179). Involved in the endosome-to-plasma membrane trafficking and recycling of SNX27-retromer-dependent cargo proteins, such as GLUT1 (PubMed:25278552). Involved in the regulation of cytokinesis; the function may involve PTPN13 and GIT1 (PubMed:23108400). Plays a role in the formation of cilia (PubMed:27767179). Involved in cargo protein localization, such as PKD2, at primary cilia (PubMed:27767179). Involved in the presentation of the tumor necrosis factor (TNF) receptor TNFRSF1A on the cell surface, and hence in the modulation of the TNF-induced apoptosis (By similarity).[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).