

Product datasheet for **SC310849**

GOLGA6 (GOLGA6A) (NM_001038640) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GOLGA6 (GOLGA6A) (NM_001038640) Human Untagged Clone
Tag:	Tag Free
Symbol:	GOLGA6
Synonyms:	GLP; GOLGA6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC310849 representing NM_001038640.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTGGCCCCAACCTACCTCCCTCCCCACCCCATGATGTTAGAAGAATCTCGACAGAATAAATTGGCA
GCAGCCAAGAAAAAGCTAAAAGAATATCAGCAAAGGAAGAGCCCTGGTATTCCAGCAGGAGCAAGACA
AAAAAGAAAAAACTGACAGTAGCCCTGAGACAACCACTTCCGGTGGTTGCCACTCACCTGGGGATAGC
CAGTACCAAGAACTAGCAGTAGCCCTGGAGTCAAGCTCAGTGACAATCAGTCAACTCAATGAAAACATA
GAATCATTGAAACAGCAGAAGAAACAAGTGAACATCAGCTGGAAGAAGCAAGAAAAACAACATGAA
ATACACAAAGCACAAATGGAGCGGTTAGAGACAATCAACATCCTCACATTGGAAGGAGAGACTTGAAG
ACCACCTTTACCATACTAAACGTGCTGCCGACACTTGAAGAAGAGTCCAAGGATCTGGCTGGCCGC
CTGCAATACTCCTACAGCGTATTCAAGAATTGGAGCGGCTCTCTGTGTGTGTCTACACAGCAGCAG
GAAGAGGACAGGTCCTCAGCTGCAGAGAAGCGGTCTCCAGCGGTGGTTACAGCAGACCATAAAGGAG
CGGGCGTGTCTGAACGCACAGTGACACAGGTGACAGAGTCACTAAAACAAGTCCAGCTAGAGCGAGAC
GAATATGCTAAACACATAAAAGGAGAGAGGGCCCGGTGGCAGGAGAGGATGTGGAATATGTCGGTGGAG
GCTCGAACATTGAAGGAAGAGAAGAAGCGTGACATACATCGGATACAGGAGCTGGAGAGGAGCTTGTC
GAATCAAAAACAGATGGCTGAGCCCCCATCCCTGGCGCCCCAGCAGTGACCTCTGTGGTGAACAG
CTACAAGATGAGGCCAAACACCTGAGGCAGGAGGTGGAAGGTCTGGAGGAAAGCTCCAGTCCCAGGTG
GAAAACATCAGGCCCTTGAAGTCTCTTAGCAAGGAACAAAGCAGAGACTCCAGGAGCAGGAGGAGATG
CTCCGAGAGCAGGAGGCGCAGAGAGTGCGGGAGCAGGAGAGAGTGTGTGAACAAACAGAGAGGCTTCGG
GAGCAGCAGAAGACGCTACAGGAGCAGGCTGAGAGGCTGCGAAAGCAGGAGCAGAGGCTACGCAACAG
GAGGAGAGGCTGCGAAAGGAGGAGGAGGCTGCAAAAGCAGGAAAGAGGCTGTGGACAGGAGGAG
AGGCTGTGGAAGAAGGAGGAGAGGCTACAAAAGCAGGAGGAGAGGCTCGCGCTCTCCAGAACCAAG
CTCGACAAGCAGCTGGCCGAGCCACAGTGCAGCTTCGAGGATCTGAACAACGAGAAAAAGAGCGACTG
CAGTTGGAGCAGCAAGTAAAGGAGCTGCAGGAGAAGCTAGACGAGGAGCACCTAGAAGCTGCCAGCCAC
CAGAACCAACAGCTAGAGACCCAGCTAAGCCTCGTGGCTCTCCCTGGAGAAGGAGATGGAGGACAACAT
CTGGACAGTGAGGAGGAGGAGGCGCTCGGCCACGCCAAACATCCCAGAGGACCTGGAGAGCCGGGAG
GCCACGAGCAGCTTTATGGACCTCCGAAGGAGAAGGCGGACGGGACGAGCAGGTGGAGAGACGAGAG
CTTGATTCTGTCAGCCTTCTGGAGTGACAGACGGCATGAGAGAGTCTTACCGTATATGAAAGCCAG
GGGGCAGTGCCAAACACGCGGCACCAGGAGATGGAGGATGTCATCAGGCTGGCCAGAAGGAGGAGGAG
ATGAAGGTGAAGCTGCTGGAGCTGCAAGAGTTGGTGTGCCCCCTGTGGGCAACCATGAGGGGCATGGC
AAATTCCTCATCGCTGCCAGAACCCTGCTGATGAGCCCACTCCAGGGGCCCCAGCCCCCAGGAATT
GGGGCTGCCGGTGAGCAGGATGTTTTTATGAAGTGAGCCTGGACAACAACGTGGAGCCTGCACCAGGA
GCGGCCAGGGAGGTTCTCCCATGACAACCCACTGTACAGCAGATCGTGACAGTGTCTCTGTCTCATG
CAGGACACCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001038640
Insert Size: 2082 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001038640.2</u>
RefSeq Size:	3176 bp
RefSeq ORF:	2082 bp
Locus ID:	342096
UniProt ID:	<u>Q9NYA3</u>
Cytogenetics:	15q24.1
MW:	79.9 kDa
Gene Summary:	The Golgi apparatus, which participates in glycosylation and transport of proteins and lipids in the secretory pathway, consists of a series of stacked cisternae (flattened membrane sacs). Interactions between the Golgi and microtubules are thought to be important for the reorganization of the Golgi after it fragments during mitosis. The protein encoded by this gene is a member of the golgin family of proteins, whose members localize to the Golgi. This gene is found in a large, low copy repeat sequence or duplicon that is found in multiple copies, that are greather than 90% similar, on chromosome 15. Duplicons are associated with deletions, inversions and other chromosome rearrangements that underlie genomic disease. The protein encoded by this gene is thought to be a functional golgin protein while the majority of the related copies of this gene are thought to be transcribed pseudogenes. [provided by RefSeq, Jul 2008]