

Product datasheet for **SC310860**

ZNF718 (NM_001039127) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	ZNF718 (NM_001039127) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF718
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC310860 representing NM_001039127. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAACCTTAAACATTCAAGGATGTGGCCATAGAATTCTCTCCAGAAGAGTGGAAATGTCTGGACACT
TCCCAGCAGAATTTATATAGAGATGTGATGTTGGAGAAGTACAGAAACCTGGTCTCCCTGGGTGTTAGT
ATCTCTAACCCAGACCTGGTCACCAGTCTGGAGCAAAGAAAGAGCCCTACAATTTGAAGATACATGAA
ACAGCAGCCAGACCCCGAGCTGTGTGTTCTCATTTACCCAAAACCTTTGGACAGTGCAGGGCATAGAA
GATTCATTCCACAACTTATACAAAAGGACATGAGAAACGTGGACATGAGAATTTAAGAAAACTTGT
AAAAGTATAAATGAGTGTAAAGTGCAGAAAGTGGTTATAATAGAATTAACCAATGCTTATTAACCTACC
CAGAAAAAAACAATTCAATCTAATATATGTGCAAGTTTTTCATAAATTTCAAATTTCAAACAAAGAT
AAGATAAGATATACTGGAGATAAAACCTTTAAATGTAAGAATGTGGCAAATCATTTCACGTGCTCTCA
CGCCTAACTCAACACAAAAGAATTCATACTGGAGAGAACCCTACACATGTGAAGAATGTGGCAAAGCC
TTAATTGGTCCCTCAATTCTTAAACATAAGAGAATTCATGCCAGAGAGAAAATTCACAAGTGTGAA
GAATGTGGTAAAGGCTTTACTCGGTCTCACACCTTAAACATAAGAGAATTCATACTGGAGAGAAA
CCCTACATATGTGAAAAATGTGGTAAAGCTTTAAACCAATCCTCAACCCTAATTTACATAAGAGAATT
CATTCTGCACAAAATACTACAAATGTGAAGAATGTGGTAAAGCCTTTAAGTGGTCTCATCCCTAAT
GAACATAAGAGAATTCATGCTGGAGAGAAACCCTTCTCATGCGAAGAATGTGGCAATGTCTTTACCACA
TCCTCAGACTTTGCTAAACATAAGAGAATTCATACAGGAGAGAAACCCTACAAATGTGAAGAATGTGGA
AAATCCTTTAATAGGTCCACAACCTTACGACACATAAGAGAATCCATACTGGAGAGAAACCCTACACA
TGTGAAGAATGTGAAAAAGCCTTTAATTGGTCTCAACCCTAATGTACACAAGAGAATTCCTCTGGA
AAAAATCCCTACAAATGTGAAGATTGTGGCAAAGCCTTTAAAGTGTGGCAAACCTGCATAATCATAAG
AAAATTCATACTGGAGAGAAACCCTACATATGTAACAATGTGGCAAAGCCTTTAAACAGTCTCACAC
TTGAATAAACATAAGAAAATTCACACTGTAGATAAACCCTACAAATGTAAAGAATGCGGAAAGCTTTT
AAGCAGTACTCCAACCTTCTCAACATAAGAGAATCATACTGGAGAAAATTTAG
ACGGGTACGGGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001039127

Insert Size: 1437 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001039127.5](#)

RefSeq Size: 3669 bp

RefSeq ORF: 1437 bp

Locus ID: 255403

UniProt ID: [Q3SXZ3](#)

Cytogenetics: 4p16.3

MW: 55.4 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) encodes the longer isoform (a).