

Product datasheet for **SC310874**

KIAA0317 (AREL1) (NM_001039479) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KIAA0317 (AREL1) (NM_001039479) Human Untagged Clone
Tag:	Tag Free
Symbol:	AREL1
Synonyms:	FIEL1; KIAA0317
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC310874 representing NM_001039479.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGTTTTACGTTATTGGTGAATCACAGTGTCTGTGGTTGCATTCTTCTTCAACAATTAAGTTCCTCTTT
GAGCTTGCCGCACGTGTAGTCAGCTTCCCTCCAGAATGAGGACCGGAGCGCCGAGGGACCGGACTATT
TATGACTACGTGCGGGGAAATTACCTGGATCCCGGTCTTGCAAAGTCTCCTGGGATTGGAAGACCCC
TATGAGGTGGGCCACAGCATGGCCTTCCGAGTGCATTTATTCTATAAGAACGGGCAGCCTTTCCTGCA
CATCGGCTGTGGGACTAAGAGTTCACATCTCTCATGTGAGCTAGCAGTGGAAATTCAGTGACCCAG
GAAGTCTTCAGGAGCCCAATCCAACGTAGTAAAAGTGGCCTTCACTGTGCGCAAGGCTGGGCGTTAT
GAAATCACAGTGAAGCTTGGTGGATTAATGTGGCATATAGTCCCTACTACAAAATTTTTCAACCTGGA
ATGGTGGTTCTTCTAAGACAAAATTGTGTGCCACTTTTCTACTCTTGATTGACCTGTGGGACGCCG
CACACCTTCAAATAGTACCCGAGATGAGTATGATAATCCACCAACAATTCATGTCTTGAGAGAT
GAGCACAATTACACCTTGTCCATTCATGAGCTCGGCCCTCAAGAAGAAGAGAGTACTGGTGTCTCATTT
GAGAAATCAGTAACATCCAACAGGCAGACTTTCAGGTGTTCTTGCGACTCACCTGCATTCTCGAGGC
TGCTTCCATGCTTGCATTTACACAAAATCAGCCAATCAATAATGGTGAATTTGACATTATTGTCTA
AGTGAGGATGAGAAGAATATCGTCGAACGCAATGTGTCCACTTCAGGCGTGAGCATTACTTTGAGGCT
TATCTTTATAATGCTACCAACTGTAGCAGCACTCCATGGCACCTGCCACCCATGCACATGACCTCTCC
CAGCGCCGGCCATCCACTGCTGTTGACGAGGAAGATGAAGACTCGCCCTCTGAGTGCCACACCCCTGAG
AAGGTGAAGAAACCGAAGAAGGTGTACTGCTATGTGTACCAAAGCAATTCAGTGAAGGAGTTCTAC
CTGAAGATCATCCCTGGCGCTTTACACCTCCGAGTGTGTCCAGGAACAAAATTTTATACCTTGGT
CAGTACCTGTCCATAAGCTGTCACACTGGTGGTGGATGATGGCATTCAACCTCTGTGGAGCTCAGC
TGTAAGGAGAGGAACATTCTAGCAGCCACTTTTATCCGCTCCCTGCATAAGAACATAGGAGGCTCTGAG
ACCTTTCAGGACAAGGTGAACTTTTTCCAGCGAGAGCTTCGCGAGGTACATATGAAAAGACCACATTCC
AAAGTCACCCTGAAGGTGAGCAGACATGCCTTGTGGAATCGTCTCTGAAAGCCACTCGGAATTTCTCC
ATCTCAGATTGGAGCAAGAACTTTGAGGTTGTTTTCCAGGATGAAGAAGCTCTGGACTGGGGAGGGCCT
CGCCGGGAATGGTTTGGCTAATCTGCAAAGCACTATTTGATACCACCAATCAGCTCTTACCCGGTTC
AGTGACAACAACCAAGCATTAGTGCATCCCAACCCTAATCGCCCCGCTCATCTGCGCCTGAAAATGTAT
GAGTTTGGGGACGGCTCGTGGCAAGTGTCTCTATGAGTCTCTTAGGAGGAGCTACAAGCAGTTG
GTCCGAGCTCGTTCACCCGCTTTTCTGGCCCAATCATAGGACTGCGTATGCATTACAAGTACTTT
GAAACAGATGACCCAGAATTCTACAAATCTAAAGTTGTTTTATCCTCAACAATGACATGAGTGAGATG
GAGCTGGTCTTTGCAGAAGAGAAATATAATAAATCAGGTCAATTGGATAAGGTTGTAGAACTCATGACA
GGTGGAGCTCAAACCTCCAGTCACCAATGCGAATAAAATCTTCTATTTAAATTTGCTGGCCCAATATCGG
CTGGCCAGTCAAGTGAAGAGGAGGTGGAACATTTCTAAAAGGCTGAATGAATTGGTCCCTGAGAAC
CTTTTGGCTATTTTTGATGAGAAATGAGCTTGTGCTGATGTGTGGGACTGGAGACATCAGTGTGTCT
GACTTCAAAGCCCATGCAGTAGTTGTTGGTGGCTCATGGCATTTCAGAGAAAAGGTCATGAGGTGGTTT
TGGACTGTGGTTTCCAGTCTGACCCAGGAGGAGTTGGCTCGGCTACTTCAGTTCAACAAGGCTCCTCT
CAGCTACCACCTGGAGGCTTTGCCGCCCTGTCCCTCATTTTCAAGATTATTGCCGCTCCGACCCATAGC
ACGCTGCCTACTGCACACACATGTTTTAACAGCTGTGCCTCCCTACATATGACTCCTATGAAGAGGTG
CACAGGATGCTGCAGCTGGCCATCAGCGAGGGTTGCGAGGGCTTTGGCATGCTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001039479

Insert Size: 2472 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:	<u>NM_001039479.1</u>
RefSeq Size:	5445 bp
RefSeq ORF:	2472 bp
Locus ID:	9870
UniProt ID:	<u>O15033</u>
Cytogenetics:	14q24.3
Protein Families:	Druggable Genome, Transmembrane
MW:	94.2 kDa
Gene Summary:	E3 ubiquitin-protein ligase which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates. Inhibits apoptosis by ubiquitinating and targeting for degradation a number of proapoptotic proteins including DIABLO/SMAC, HTRA2 and SEPT4/ARTS which are released from the mitochondrion into the cytosol following apoptotic stimulation.[UniProtKB/Swiss-Prot Function]