

Product datasheet for **SC310898**

ATP6V1C2 (NM_001039362) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ATP6V1C2 (NM_001039362) Human Untagged Clone
Tag:	Tag Free
Symbol:	ATP6V1C2
Synonyms:	ATP6C2; VMA5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC310898 representing NM_001039362. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCCGGAGTTTTGGTTAATTTCTGCCCTGGCGATAAGGAAAATTTGCAAGCTCTGGAGAGGATGAAT
ACTGTAACCTCCAAGTCCAACCTGTCTTATAATACCAAATTCGCTATTCTGACTTCAAGGTGGGGACC
TTGGATTCCCTGGTTGGCTCTCTGATGAGTTGGGGAACTCGACACCTTTGCTGAAAGCCTCATAAGG
AGAATGGCTCAGAGCGTGGTGAAGTCATGGAGGACTCAAAGGGGAAGGTCCAGGAGCACCTCCTGGCA
AACGGAGTTGACTTAACATCCTTTGTGACCCACTTTGAATGGGACATGGCCAAATATCCTGTCAAGCAG
CCGCTCGTGAGTGTGGTGGACACAATAGCCAAGCAACTGGCGCAGATCGAGATGGACCTGAAGTCCCGA
ACGGCCGCTACAACACTCTGAAGACAACTGGAGAACCTGGAAAAGAAATCCATGGGGAACTCTTC
ACCCGGACACTGAGTGATATTGTGAGCAAAGAGGACTTCGTGCTGGATTCTGAATATCTCGTCAACTT
CTGGTCATCGTCCCAAACCAAACTACTCACAATGGCAAAAAACCTACGAATCTCTCAGACATGGTG
GTCCCTCGATCAACCAAACTCATTACTGAGGACAAGGAAGGGGGCCTTTTCACTGTGACTCTGTTTCGA
AAAGTGATTGAAGATTTCAAACCAAGGCCAAAGAAAACAAGTTCCTGTTCTGTAATTTTACTATGAT
GAGAAGGAAATTTGAAAGGAAAGGGAGGAGATGGCCAGATTGCTGTCTGATAAGAAGCAACAGTACAA
ACTTCTGTGTTGCTCTTAAAAAGGGATCATCCACCTTCCCGACCACAAGGTTAAGGTAACCCCGCTA
GGTAACCTGATAGGCCGTGCTGCGGGGACAGCCGACAGAGAGAGAGTGGGGCGAGGGTGAGGGC
CCCTGCTGCGCTGGCTCAAGGTGAACCTCAGTGAAGCCTTCATTGCTGGATCCACATCAAGGCCCTG
AGAGTGTGTTGAGTCCGTGCTCAGGTATGGACTACCAAGTGAACCTCCAGGCAGTGCTCCTGCAGCCG
CATAAGAAGTCATCCACCAAGCGTTTAAGAGAGGTTCTAAACTCTGTCTCCGACATCTGGATGAAGTA
GCCGCTACAAGTACTGGATGCATCTGTGGAGATCCCGGGACTGCAACTCAATAACCAAGACTATTTT
CCTTATGTCTACTTCCATATTGACCTTAGTCTTCTTGACTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001039362

Insert Size: 1284 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001039362.1](#)

RefSeq Size: 3250 bp

RefSeq ORF: 1284 bp

Locus ID: 245973

UniProt ID:	<u>Q8NEY4</u>
Cytogenetics:	2p25.1
Protein Pathways:	Epithelial cell signaling in Helicobacter pylori infection, Metabolic pathways, Oxidative phosphorylation, Vibrio cholerae infection
MW:	48.8 kDa
Gene Summary:	<p>This gene encodes a component of vacuolar ATPase (V-ATPase), a multisubunit enzyme that mediates acidification of eukaryotic intracellular organelles. V-ATPase dependent organelle acidification is necessary for such intracellular processes as protein sorting, zymogen activation, receptor-mediated endocytosis, and synaptic vesicle proton gradient generation. V-ATPase is composed of a cytosolic V1 domain and a transmembrane V0 domain. The V1 domain consists of three A, three B, and two G subunits, as well as a C, D, E, F, and H subunit. The V1 domain contains the ATP catalytic site. This gene encodes alternate transcriptional splice variants, encoding different V1 domain C subunit isoforms. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (a).</p>