

Product datasheet for **SC310962**

DCAMKL2 (DCLK2) (NM_001040260) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DCAMKL2 (DCLK2) (NM_001040260) Human Untagged Clone
Tag:	Tag Free
Symbol:	DCLK2
Synonyms:	CL2; CLICK-II; CLICK2; CLIK2; DCAMKL2; DCDC3; DCDC3B; DCK2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC310962 representing NM_001040260.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCCAGCACCAGGAGTATCGAGCTGGAGCACTTTGAGGAACGGGACAAAAGGCCGCGCCGGGTTCG
CGGAGAGGGGCCCCAGCTCCTCCGGGGCAGCAGCAGCTCGGGCCCAAGGGGAACGGGCTCATCCCC
AGTCCGGCGCACAGTGCCCACTGCAGCTTCTACCGCACGCGGACCCTGCAGGCCCTCAGCTCGGAGAAG
AAGGCCAAGAAGGCGCGTCTACCGGAACGGGACCGCTACTTCAAGGGCCTGGTGTTCATCTCC
AGCGACCGCTTCCGGTCTTCGATGCGCTCCTCATAGAGCTCACCCGCTCCCTGTCGGACAACGTGAAC
CTGCCCCAGGGTGTCCGCACTATCTACACCATCGACGGCAGCCGGAAGGTACCAGCCTGGACGAGCTG
CTGGAAGGTGAGAGTTACGTGTGTGCATCCAATGAACCATTTTCGTAAGTCGATTACACAAAAATATT
AATCCAACTGGTCTGTGAACATCAAGGGTGGGACATCCCAGCGCTGGTGTGCCTCCTCTGTGAAA
AGTGAAGTAAAAGAAAGTAAAGATTTCAAAACCAAGTTAGTGACTGTGATTCGAAGTGGAGTGAAG
CCTAGAAAAGCCGTGCGGATCCTTCTGAATAAAAAGACTGCTCATTCTTTGAACAAGTCTTAACAGAT
ATCACCGAAGCCATTAAGTACTAGACTCAGGAGTCGTAAGAGGCTCTGCACCCTGGATGAAAGCAGGTT
ACTTGTCTGCAAGACTTTTTTGGTGATGACGATGTTTTATTGCATGTGGACCAGAAAAATTCGTAT
GCCAAGATGACTTTGTCTGGATCATAGTGAATGTCGTGTCCTGAAGTCATCTTATTCTCGATCCTCA
GCTGTTAAGTATTCTGGATCCAAAAGCCCTGGGCCCTCTCGACGCAGCAAATCACCAGCTTCAGTTAAT
GGAAGTCCCAGCAGCCAACTTTCTACTCTAAATCTACGAAATCCTCCAGTTCCTCTCCAAGTATGCA
GGAAGTTTCAGAGGATTAAGCAGATTTCTGCTCATGGCAGATCTTCTTCCAATGTAACGGTGGACCT
GAGCTTGACCGTTGCATAAGTCTGAAGGTGTGAATGGAACAGATGCTCTGAATCATCAACTCTTCTT
GAGAAATACAAAATTGGAAGGTCAATGGTGATGGCAATTTTTCAGTAGTCAAAGAGTGTATAGCAGG
TCCACTGGAAAGGAGTTTGCCCTAAAGATTATAGACAAAAGCCAAATGTTGTGAAAGGAACACCTGATT
GAGAATGAAGTGTCAATACTGCGCCGAGTGAACATCCCAATATCATTATGCTGGTTCGAGGAGATGGAA
ACAGCAACTGAGCTCTTTCTGGTGATGGAATTGGTCAAAGGTGGAGATCTCTTTGATGCAATTAATTCTG
TCGACCAAGTACACTGAGAGAGATGGCAGTGCCATGGTGTACAACCTAGCCAATGCCCTCAGGTATCTC
CATGGCCTCAGCATCGTGCACAGAGACATCAAACCAGAGAATCTCTTGGTGTGTGAATATCCTGATGGA
ACCAAGTCTTTGAAACTGGGAGACTTTGGGCTTGGCAGTGTGGTGAAGGCCCTTATACACAGTCTGT
GGCACACCCACTTATGTGGCTCCAGAAATCATTGCTGAAACTGGCTATGGCTGAAGGTGGACATTTGG
GCAGCTGGTGTGATCACATACATACTTCTCTGTGGATTCCCACCATTCCGAAGTGAACAATCTCCAG
GAAGATCTCTCGACCAGATCTTGGCTGGGAAGCTGGAGTTTCCGGCCCCCTACTGGGATAACATCACG
GACTCTGCCAAGGAATTAATCAGTCAAATGCTCAGGTAATGTTGAAGCTCGGTGTACCGCGGGACAA
ATCCTGAGTACCCCTGGGTGTGAGATGATGCCTCCCAGGAGAATAACATGCAAGCTGAGGTGACAGGT
AACTAAAACAGCACTTTAATAATGCGCTCCCCAACAGAACAGCACTACCACCGGGTCTCCGTATC
ATGAACACGGCTTAGATAAGGAGGGGCAGATTTTCTGCAGCAAGCACTGTCAAGACAGCGGCAGGCCT
GGGATGGAGCCCATCTCTCCAGTTCCTCCCTCAGTGGAGGAGATCCCTGTGCCTGGGGAAGCAGTCCC
GCCCCACCCCTCCGGAATCTCCCACCCCCACCCTCCTCCCGCTGCCCCGGGTGGTGGAGCGGCAGGA
ACCTGGCGCCGCCACCGAGACTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001040260

Insert Size: 2301 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001040260.3](#)

RefSeq Size: 4274 bp

RefSeq ORF: 2301 bp

Locus ID: 166614

UniProt ID: [Q8N568](#)

Cytogenetics: 4q31.23-q31.3

Protein Families: Druggable Genome, Protein Kinase

MW: 83.6 kDa

Gene Summary: This gene encodes a member of the protein kinase superfamily and the doublecortin family. The protein encoded by this gene contains two N-terminal doublecortin domains, which bind microtubules and regulate microtubule polymerization, a C-terminal serine/threonine protein kinase domain, which shows substantial homology to Ca²⁺/calmodulin-dependent protein kinase, and a serine/proline-rich domain in between the doublecortin and the protein kinase domains, which mediates multiple protein-protein interactions. The microtubule-polymerizing activity of the encoded protein is independent of its protein kinase activity. Mouse studies show that the DCX gene, another family member, and this gene share function in the establishment of hippocampal organization and that their absence results in a severe epileptic phenotype and lethality, as described in human patients with lissencephaly. Multiple alternatively spliced transcript variants have been identified. [provided by RefSeq, Sep 2010]
Transcript Variant: This variant (1) is known as exon-5-lacking transcript (PMIDs:15611072 and 18075264). It lacks an in-frame coding exon in the 5' region, as compared to variant 2. The resulting isoform (a) lacks an internal segment, as compared to isoform b. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.