

## Product datasheet for **SC310966**

### AGTRAP (NM\_001040195) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	AGTRAP (NM_001040195) Human Untagged Clone
Tag:	Tag Free
Symbol:	AGTRAP
Synonyms:	ATRAP
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001040195, the custom clone sequence may differ by one or more nucleotides ATGGAGCTGCCTGCTGTGAACCTGAAGGGGCTGCATTGTATTCTCAGGCTCCTATGCCTG GGCCAACCTCACCATCCTGGCCTTGGGCGTGTGGGCTGTGGCTCAGCGGGACTCCATCGA CGCCATAAGCATGTTTCTGGGTGGCTTGTGGCCACCATCTTCTGGACATCGTGCACAT CAGCATCTTCTACCCGCGGGTCAGCCTCACGGACACGGGCGCTTGGCGTGGGCATGGC CATCCTCAGTTGCTGCTCAAGCCGCTCTCCTGCTGCTTCGTCTACCACATGTACCGGA GCGCGGGGTTTCCTTGGGTCTTCTCAGGACCGTAG
Restriction Sites:	Please inquire
ACCN:	NM_001040195
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001040195.1](#), [NP\\_001035285.1](#)

**RefSeq Size:** 1144 bp

**RefSeq ORF:** 336 bp

**Locus ID:** 57085

**UniProt ID:** [Q6RW13](#)

**Cytogenetics:** 1p36.22

**Protein Families:** Druggable Genome, Transmembrane

**Gene Summary:** This gene encodes a transmembrane protein localized to the plasma membrane and perinuclear vesicular structures. The gene product interacts with the angiotensin II type I receptor and negatively regulates angiotensin II signaling. Alternative splicing of this gene generates multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]

Transcript Variant: This variant (3) lacks an exon and uses an alternate splice site for a second exon in the 3' coding region, compared to variant 1, resulting in a frameshift. The resulting protein (isoform c) has a distinct C-terminus compared to isoform a.