

Product datasheet for **SC310982**

TBC1D3G (NM_001040282) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TBC1D3G (NM_001040282) Human Untagged Clone
Tag:	Tag Free
Symbol:	TBC1D3G
Synonyms:	TBC1 domain family, member 3G
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC310982 representing NM_001040282.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGACGTGGTAGAGGTCGCGGGTAGTTGGTGGGCACAAGAGCGAGAGGACATCATTATGAAATACGAA
AAGGGACACCGAGCTGGGCTGCCAGAGGACAAGGGGCCTAAGCCTTTTCGAAGCTACAACAACAACGTC
GATCATTGGGGATTGTACATGAGACGGAGCTGCCTCCTCTGACTGCGCGGGAGGCGAAGCAAATTCGG
CGGGAGATCAGCCGAAAGAGCAAGTGGTGGATATGCTGGGAGACTGGGAGAAATACAAAAGCAGCAGA
AAGCTCATAGATCGAGCGTACAAGGGAATGCCATGAACATCCGGGGCCCGATGTGGTCAGTCCTCTG
AACATTGAGGAAATGAAGTTGAAAAACCCGGAAGATACCAGATCATGAAGGAGAAGGGCAAGAGGTCA
TCTGAGCACATCCAGCGCATCGACCGGGACATAAGCGGGACATTAAGGAAGCATATGTTCTTCAGGGAT
CGATACGGAACCAAGCAGCGGAACTACTCCACATCCTCTGGCATATGAGGAGTATAACCCGGAGGTG
GGCTACTGCAGGGACCTGAGCCACATCGCCGCCTTGTCTCTCTATCTTCTGAGGAGGATGCATTC
TGGGCACTGGTGCAGCTGCTGGCCAGTGAGAGGCACTCCCTGCAGGGATTTCACAGCCCAAATGGCGGG
ACCGTCCAGGGGCTCCAAGACCAACAGGAGCATGTGGTAGCCACGTACAAACCAAGACCATGGGGCAT
CAGGACAAGAAAGATCTATGTGGGCACTGTTCCCGTTAGGCTGCCTCATCCGGATATTGATTGACGGG
ATCTCTCTCGGGCTCACCTGCGCTGTGGGACGTGTATCTGGTAGAAGGCGAACAGGCGTTGATGCCG
ATAACAAGAAATCGCCTTTAAGGTTACAGAGAAGCGCCTCACGAAGACGTCCAGGTGTGGCCCGTGGGCA
CGTTTTTGCAACCGGTTCTGTGATACCTGGGCCAGGGATGAGGACACTGTGCTCAAGCATCTTAGGGCC
TCTATGAAGAACTAACAAAGAAAGCAGGGGACCTGCCACCCAGCCAAACCCGAGCAAGGGTCGTCG
GCATCCAGGCCTGTGCCGGCTTACGTGGCGGGAAGACCCTCTGCAAGGGGACAGGCAGGCCCTCCA
GGCCACCAAGCCCGTTCCCGGCCCATTTGGTCAGTTCCCGGCCACGGGCACCTCGTTCTTCCACA
CCCTGTCTGTGGGGCTGTCCGGGAAGACACCTACCCTGTGGGCACTCAGGGTGTGCCAGCCCGGCC
CTGGCTCAGGGAGGACCTCAGGGTTCCTGGAGATTCTGCAGTGGAACCTCATGCCCGCCTCCCAACG
GACCTGGACGTAGAGGGGCTTGGTTCCGCCATTATGATTTAGACAGAGCTGCTGGGTCCGTGCCATA
TCCAGGAGGACAGCTGGCCCCCTGCTGGCAGGCTGAACACCCTGCGGAGCGGGTGAGATCGGCTTTC
GCTGCACCCAGCACTGATTCCGACAGGGCACCCCTTCAGAGCTAGGGACGAACAGCAGTGTGCTCC
ACCTCAGGGCCTTGCTCTGCGGCTCCACTTGAAAGTTCTCAGTTCCCTCCAGGCTTCAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001040282

Insert Size: 1650 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001040282.2](#)

RefSeq Size: 2177 bp

RefSeq ORF: 1650 bp

Locus ID: 654341

Cytogenetics: 17q12

MW: 62.2 kDa