

Product datasheet for SC310999

Lunatic Fringe (LFNG) (NM_001040168) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	Lunatic Fringe (LFNG) (NM_001040168) Human Untagged Clone
Tag:	Tag Free
Symbol:	LFNG
Synonyms:	SCDO3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC310999 representing NM_001040168. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTCAAGCGCTGCGGCCGCGCCTGCTGCTGGCGCTGGCGGGCGCGCTGCTCGCCTGCCTGCTGGTG
CTCACCGCCGACCCGCGCCGCCCTCCACTGCCCGCGAGCGCGCCGGCGCGCTGCGCAGCCTGGCG
GGCCCCGCGGGGGCTGCCCGCGCCCGGGCTGGGGGCGCGCGCGCGCCCGGGGGCGCTGGTCCGC
GACGTGCACAGTCTGTCCGAGTACTTCAGCCTGCTCACCCGCGCGCGCAGAGATGCGGGCCCGCCGCC
GGGGTGGCCCCCGCCCGACGGCCACCCGCGCCCTGGCCGAGCCGCTCGCGCCCGAGACGTC
TTCATCGCTGTCAAGACCACAAAAGTTCACCGCGCGCGCCTCGACCTGCTGTGGAGACCTGGATC
TCGCGCCACAAGGAGATGACGTTTCATCTTCACTGACGGGGAAGATGAGGCCCTGGCCAGGCACACGGC
AACGTGGTCATCAAACTGCTCGGCCGCCACAGCCGCCAGGCGCTGTCTGCAAGATGGCCGTGGAG
TATGACCGCTTACGAGTCCGGCAGGAAGTGGTCTGCCACGTGGACGATGACAACCTACGTCACCTG
CGGGCCCTGCTGCGGCTGCTGGCCAGCTACCCGCACACGCGGGACGTCTACGTGCGCAAGCCAGCCTG
GACAGGCCATCCAGGCCATGGAGCGGGTACGCGAGAACAAGGTGCGTCTGCTCCACTTCTGTTTGC
ACGGGCGCGCTGGCTTCTGCATCAGCCGTGGGCTGGCTCTGAAGATGAGCCCGTGGCCAGCGGGGT
CACTTCATGAATACGGCTGAGCGGATCCGGCTGCCTGATGACTGCACCATCGGCTACATCGTGGAGGCC
CTGCTGGGTGTGCCCTCATCCGACGCGCCTTCCACTCCACCTGGAGAACCCTGCAGCAGGTGCC
ACCTCGGAGTCCACGAGCAGGTGACGCTGAGCTACGGTATGTTTGAAAACAAGCGGAACGCCGTCAC
GTGAAGGGGCCCTTCTCGGTGGAGGCCGACCCATCCAGGTGGGGAACCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001040168

Insert Size: 1086 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001040168.1](#)

RefSeq Size: 1403 bp

RefSeq ORF: 1086 bp

Locus ID: 3955

UniProt ID:	Q8NES3
Cytogenetics:	7p22.3
Protein Families:	Transmembrane
Protein Pathways:	Notch signaling pathway
MW:	39.6 kDa
Gene Summary:	<p>This gene is a member of the glycosyltransferase 31 gene family. Members of this gene family, which also includes the MFNG (GeneID: 4242) and RFNG (GeneID: 5986) genes, encode evolutionarily conserved glycosyltransferases that act in the Notch signaling pathway to define boundaries during embryonic development. While their genomic structure is distinct from other glycosyltransferases, these proteins have a fucose-specific beta-1,3-N-acetylglucosaminyltransferase activity that leads to elongation of O-linked fucose residues on Notch, which alters Notch signaling. The protein encoded by this gene is predicted to be a single-pass type II Golgi membrane protein but it may also be secreted and proteolytically processed like the related proteins in mouse and Drosophila (PMID: 9187150). Mutations in this gene have been associated with autosomal recessive spondylocostal dysostosis 3. [provided by RefSeq, May 2018]</p> <p>Transcript Variant: This variant (2) includes an alternate 3' terminal exon compared to variant 1. The resulting protein (isoform b) is shorter and has a distinct C-terminus compared to isoform 1.</p>