

## Product datasheet for **SC311017**

### Striatin 4 (STRN4) (NM\_001039877) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Striatin 4 (STRN4) (NM_001039877) Human Untagged Clone
Tag:	Tag Free
Symbol:	Striatin 4
Synonyms:	PPP2R6C; ZIN; zinedin
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC311017 representing NM\_001039877.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGATGGAGGAGCGAGCGGCCCGCCGGTCGCCCGCCCGCTCCTCCTGCCGTCCGCTCGGCTCAGGC
GGGGCCCTGGCCCACCTGGGGCGGCCCGGTCTCCGCCCTGCCCGGGCCGGCCCGCAGGTAAG
GGAGGCGGCGGGAGGAGCCCGGGCCACGGCGGGCCGGAGCCCTGAGCCTGCCGGGATCCTG
CACTTTATCCAGCACGAGTGGGCGCTTCAAGCCGAGAAAGCCGCTGGGAGGCCGAGCGCGCCGAG
TTACAGGCTCAGGTGGCCTTCTTCAGGGAGAGAGAAAGGGCAGGAGAATCTAAAGACGGACCTGGTG
CGGCGGATCAAGATGCTAGAGTATGCGCTGAAGCAGGAAAGGGCCAAATATCATAAACTGAAGTTTGGG
ACAGACCTGAACAGGGGAGAAGAAAGCAGATGTGTGAGAACAAGTCTCAATGGCCCCGTGGAATCG
GTCACCCTGGAGAACAGCCGTTGGTGTGAAGGAGGGGGCAGCTTCTCCGACGTACCTGGAAGAG
GTGGGCTACACAGACACCATCTCGACATGCGGTCCAAGCGCTCCGTTCCCTGCTGGGCCGCTCGCTG
GAGCTCAACGGGCGAGTGGAGCCGAGTGAAGGGCCCCAGGGTCCACCAGGCCCTGCAGGGCTCAGT
GGTGGGGAGTCCGTGCTGGTGAACAGATCGAGGAGCAGATAAAGAGGAACGCGGAGGCAAGATGGC
AAAGAGCGCTTGGGCGGCTCAGTGTGGGCGAGATCCCCTTCTGCAGAACTGCGAGGACGAAGACAGC
GACGAGGACGATGAGCTGGACAGCGTGCAGCACAAAGCAGCGTGTGAAGTCCCATCCAAGGCTCTG
GTGCCCGAAATGGAAGACGAGGATGAGGAAGACGACTCTGAGGATGCTATCAATGAGTTTGATTTCTG
GGCTCAGGAGAGGATGGGAAGGGGCTCCAGACCCTCGGCGGTGCACTGTGGATGGGAGCCCCATGAG
CTGAAAGCCGTGGGTCAAACCTCAAGGCATTCTGGTGAACCTGCGGGATGTGGATGGGCTGCCCCA
AAAGTGAAGCCCGCTCCTGGCACACCCAGCCCGGCCACATGAAGTTCTTTGGCTTCTCCTCA
GACCTTTCATCATGGACACTATCGGGGGCGGGAGGTGAGCCTGGGGACTTGGCAGATCACCAGTC
ACCAACGACAACGACCTCAGCTGCGATCTGTCTGACAGCAAAGATGCTTTTAAGAAGACGTGGAACCC
AAGTTCACCCTGCGCTCGCACTACGACGCACTTCCCTGGCCTTCCACCACAGCCAGTCCGGCTCTG
CTCACCCTCCGAGGACGGCACGCTCAAGCTCTGGAACCTGCAGAAGCGGTACCGCCAAAGAAGAT
GCGGGCTAGATGTGGAACCTATACATGCTTCCGGGCTCACAGGGGCCAGTGTGGCTGTGGCTATG
GGCAGCAACAGTGAATACTGCTACAGTGGCGGGCAGATGCTGCATCCATAGTTGGAAGATTCCAGAC
CTCAGCATGGATCCCTATGATGGCTACGACCAAGCGTGTGAGCCAGTCTGGAGGGCCACGGGGAC
GCCGTGTGGGCTGGCCTTCACTCCACCTCCAGCGCCTGGCCTCCTGTTCTGCTGATGGCACCGTC
CGCATCTGGGACCCAGCAGCAGCAGCCGGCCTGCCTCTGCACCTTCCCACAGCCAGCAGCAGGGG
GTCCCCACCTCAGTGGCCTTACCAGCACCGAGCCTGCCACATCGTGGCCTCCTCCGCTTGGCGAC
ACCGTCTTGATGACATGGAGGTTGGCAGTGCCTCCTCACGCTGGAGTCCCGGGGACGACGCGTCCA
ACCCAGATCAACCAAGTGGTGAATCATCAAACCAGCCTCTCACCATCACCGCCACGACGACAGGGGC
ATCCGCTTCTGGACAATCGGACAGGTAAGCCGGTGCCTCCATGGTTGCACACCTGGACGCACTACC
TGCTAGCCGTGGACCCCAACGGCGCATTCCTGATGTCAGGAAGCCATGACTGCTCCCTGCGTCTCTGG
AGCCTGGACAACAAAACGTGCGTGCAGGAGATCACGGCCACCGCAAGAAGCAGCAGGAGGCCATCCAC
GCTGTTGCTGCCACCCAGCAAGGCCCTCATTGCCAGTGTGGCGCTGATGCCCTGGCCAAGTCTTC
GTATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_001039877

**Insert Size:** 2283 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001039877.1</a></u>
<b>RefSeq Size:</b>	3242 bp
<b>RefSeq ORF:</b>	2283 bp
<b>Locus ID:</b>	29888
<b>UniProt ID:</b>	<u><a href="#">Q9NRL3</a></u>
<b>Cytogenetics:</b>	19q13.32
<b>MW:</b>	81.3 kDa
<b>Gene Summary:</b>	<p>Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling protein.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) represents the longer transcript and encodes the longer isoform (2).</p>