

## Product datasheet for **SC311095**

### PTPN5 (NM\_001039970) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PTPN5 (NM_001039970) Human Untagged Clone
Tag:	Tag Free
Symbol:	PTPN5
Synonyms:	PTPSTEP; STEP; STEP61
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC311095 representing NM\_001039970.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAATTATGAGGGAGCCAGGAGTGAGAGAGAGAACCACGCTGCTGATGACTCCGAGGGAGGGGCCCTG
GACATGTGCTGCAGTGAGAGGCTACCGGGTCTCCCCAGCCGATAGTGATGGAGGCACTGGACGAGGCT
GAAGGGCTCCAGGACTCACAGAGAGAGATGCCGCCACCCCTCCTCCCTCGCCGCCCTCAGATCCAGCT
CAGAAGCCACCACCTCGAGGCGCTGGGAGCCACTCCCTACTGTGAGGAGCAGCCTGTGCCTGTTGCTG
GCCTCACAGTTCTGTGCTGGAACCCACGGCCTGGCTTACTCTGGGACGTGGGGAGTCCCCAGTCTGCTG
CTGGTCTTTCTGTCCGTGGGCTGGTCTCTGTTACCACCCTGGTGTGGCACCTCTGAGGACACCCCCA
GAGCCACCACCCACTGCCCTGAGGACAGGCGCCAGTCAGTGAGCCGCCAGCCCTCCTTCACCTAC
TCAGATGGATGGAGGAGAAGATCGAGGATGACTTCTGGACCTCGACCCGGTCCCCGAGACTCCTGTG
TTTGATTGTGTGATGGACATCAAGCCTGAGGCTGACCCACCTCACTCACGTCAGTCCATGGGCTG
CAGGAGAGGAGGGTTCCAATGTCTCCCTGACCCTGGACATGTGCACTCCGGGCTGCAACGAGGAGGGC
TTTGGCTATCTCATGTCCCACGTGAGGAGTCCGCCCGGAGTACTGCTCAGCGCTCCCGTGTCCCTC
CAAGCAGAAGAGCTTCATGAAAAGGCCCTGGACCCCTTCTGCTGCAGGCGGAATTCCTTTGAAATCCCC
ATGAACCTTTGTGGATCCGAAAGAGTACGACATCCCTGGGCTGGTGCAGGAAACCAGTACAAAACATA
CTTCCCAACCCTCACAGCAGAGTGTGTCTGACCTCACCAGACCCTGACGACCCTCTGAGTTCCTACATC
AATGCCAACTACATCCGGGGCTATGGTGGGGAGGAGAAGGTGTACATCGCCACTCAGGGACCCATCGTC
AGCACGGTCGCCGACTTCTGGCGCATGGTGTGGCAGGAGCACACGCCATCATTGTGATGATCACCAC
ATCGAGGAGATGAACGAGAAATGCACCGGAGTATTGGCCGGAGGAGCAGGTGGCGTACGACGGTGTGAG
ATCACTGTGCAGAAAGTCATTACACGGAGGATTACCGGCTGCGACTCATCTCCCTCAAGAGTGGGACT
GAGGAGCGAGGCCTGAAGCATTACTGGTTCACATCCTGGCCCGACCAGAAGACCCAGACCCGGGCCCCC
CCACTCTGCACCTGGTGCAGGAGTGGAGGAGGAGCCAGCCAGCAGGAGGGGCCCCACTGTCCCCCCTC
ATCGTCCACTGCAGTGCAGGGATTGGGAGGACCGGCTGCTTCAATGCCACCAGCATCTGCTGCCAGCAG
CTGCGGCAGGAGGGTGTGGTGGACATCCTGAAGACCAGTGCCAGCTCCGTGAGGACAGGGGCGGCATG
ATCCAGACATGCGAGCAGTACCAGTTTGTGCACCACGTCATGAGCCTCTACGAAAAGCAGCTGTCCCAC
CAGTCCCCAGAATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

**Restriction Sites:**

Sgfl-MluI

**Plasmid Map:**



<b>ACCN:</b>	NM_001039970
<b>Insert Size:</b>	1602 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001039970.1</a>
<b>RefSeq Size:</b>	3105 bp
<b>RefSeq ORF:</b>	1602 bp
<b>Locus ID:</b>	84867
<b>UniProt ID:</b>	<a href="#">P54829</a>
<b>Cytogenetics:</b>	11p15.1
<b>Protein Families:</b>	Druggable Genome, Phosphatase, Transmembrane
<b>Protein Pathways:</b>	MAPK signaling pathway
<b>MW:</b>	60 kDa
<b>Gene Summary:</b>	<p>May regulate the activity of several effector molecules involved in synaptic plasticity and neuronal cell survival, including MAPKs, Src family kinases and NMDA receptors. [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (3) has an alternate splice site in the 5' coding region, compared to variant 1. The resulting isoform (b) lacks an internal segment, compared to isoform a.</p>