

## Product datasheet for **SC311147**

### DTX3L (NM\_138287) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	DTX3L (NM_138287) Human Untagged Clone
Tag:	Tag Free
Symbol:	DTX3L
Synonyms:	BBAP; RNF143
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC311147 representing NM\_138287.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGCCTCCACCTGCGCCCGCGTCCCGCTCCTCGTGGGGTGTACAAGTCCGGCCCCGAGTACGA
AGGAAGCTGGAGAGCTACTTCCAGAGCTAAGTCTCGGGCGGGGGAGTGCACGGTCAACCCAG
GAACACGAAGCCCGGGACCTTCCGGGTGGAGTTCAGTGAAAGGCAGCTAAGGAGAGAGTGTGAAA
AAAGGAGAGCACCAAACTTGTGACGAAAACTGTGCCATTTTCTGGTACCCACTGAAAATTCA
ATAAAGAAGAACACGAGACCTCAAATTTCTTCACTGACACAATCACAAAGCAGAAACACCGTCTGGTGT
ATGCATCAACATGAAGGACATATCCTAATGCTGTGGATTCTGTCTCCAAAAGATCTTTCTTACTGTA
ACAGCTGACCTGAACTGTAACCTGTTCTCAAAGAGCAGAGGGCATAACATAACCACACTGTGCCCTAGT
ATCAGAAAAATGGAAGGTACGATGGAATTGAGAAGGTGTGTGGTACTCCAAGACATTGAAAGAATA
CATCAATTTTGGAGTGAGCAGTTCCTGGAAAGTGAGCAGAAACAACATTTTCCCTTCAATGACAGAG
AGGAAGCCACTCAGTCAGCAGGAGAGGGACAGCTGCATTTCTCCTTCTGAACCAGAAACCAAGGCAGAA
CAAAAAAGCAACTATTTTGAAGTTCCTTGCCTTACTTTGAATACTTTAAATATATCTGTCTGATAAA
ATCAACTCAATAGAGAAAAGATTTGGTGTAAACATTGAAATCCAGGAGAGTTCCTCAAATATGGTCTGT
TTAGATTTACCTCAAGTCGATCAGGTGACCTGGAAGCAGCTCGTGAGTCTTTTGTAGTGAATTTAG
AAGAACACAGAACCTCTGAAGCAAGATGTGTCTTTAGCAGACAGTAAGCAGGCAAAATAAATTTAAA
CAGGAATTGAATCACCAGTTTACAAAGCTCCTTATAAAGGAGAAAAGGAGGCGAATTAACCTCCTTGGG
ACCCAAGATGACATTTCACTGCCAAACAAAAATCTCTGAAGCTTTTGTCAAGTACCTGTGAAACTA
TTTGTCTGCAATTACATGATGAATGTAATTGAGGTTGATAGTCCCACTATAAATTTTAGAACTGAA
TTACTACAGGAGATATCAGAGATCGAAAAAGGATGACATTTGCAGCAAGGTTTCTGAGAAAGGTGAG
AAAACCTGCATTCTGTTTGAATCCAAGGACAGGCAGGTAGATCTATCTGTGCATGCTTATGCAAGTTTC
ATCGATGCCTTTCAACATGCCTCATGTGAGTGTGAGAGAAAGTCTTTTACTGAAGCTTTGGGCAAG
GAGAGAAAGCACTTACATCAGACCAAGTTTGTGATGACTTTAGAAAAAGACATCCAAATGTACACTTT
GTGCTAAATCAAGAGTCAATGACTTTGACTGGTTTGCCTTACCTTGCAGGCAAGCAGTATGTT
CTAAAAGGAGGAGGAATGTCTTTCATTGGCTGGAAAGAAATGAAAGAGGGTCAAGAACCCGATGGAC
ATTGATAGCGATGATTCAAAGCAGCTTCTCCGCACTCAAGGGCTCTGTGAGTCTGAGGCCTCAGAA
CTGGACAAGAAGGAAAAGGCATCTGTGTCATCTGTATGGACACCATTAGTAACAAAAAGTGTACCA
AAGTGCAAGCATGAATTTGCGCCCTTGTATCAACAAAGCCATGTCATATAAGCCAATCTGTCCACA
TGCCAGACTTCTATGGTATTCAGAAAGGAAATCAGCCAGAGGGAAGCATGGTTTTCACTGTTTCAAGA
GACTCACTTCCAGGTTATGAGTCTTTGGCACCATTGTGATTACTTATTCTATGAAAGCAGGCATACAA
ACAGAAGAACACCCAAACCCAGGAAAGAGATACCTGGAATACAGCGAACTGCATACTTGCCTGATAAT
AAGGAAGGAAGGAAGGTTTTGAACTGCTTTATAGGGCTTTGACCAAAAGCTGATTTTTACAGTGGGG
TACTCTCGCGTATTAGGAGTCTCAGATGTCATCACTTGGAAATGATATTCACCACAAAACATCCCGGTTT
GGAGGACCAGAAATGATGGCTATCCTGATCCTTCTTACCTGAAACGTGTCAAAGAGGAGCTGAAAGCC
AAAGGAATTGAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

**Restriction Sites:** Sgfl-Mlul

**Plasmid Map:**


**ACCN:** NM\_138287

**Insert Size:** 2223 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_138287.3](#)

**RefSeq Size:** 5768 bp

**RefSeq ORF:** 2223 bp

**Locus ID:** 151636

**UniProt ID:** [Q8TDB6](#)

<b>Cytogenetics:</b>	3q21.1
<b>Domains:</b>	RING
<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	Notch signaling pathway
<b>MW:</b>	83.6 kDa
<b>Gene Summary:</b>	DTX3L functions as an E3 ubiquitin ligase (Takeyama et al., 2003 [PubMed 12670957]). [supplied by OMIM, Nov 2009]