

Product datasheet for SC311163

KIAA0895L (NM_001040715) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KIAA0895L (NM_001040715) Human Untagged Clone
Tag:	Tag Free
Symbol:	KIAA0895L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC311163 representing NM_001040715. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTGTGGAAGCTCAGGGGCTCAGGCGTATGATCAGGCACCCCCAGCCCGCTACCACTCCCCATCC
CTGCGCCATAGGCTGAAGCCCTCAGACCGAGATGGGCCACCACTGTACCCCTGGTCTCAGTCCCTGGCC
TTGCCCTGGCTCTGGCAGTCCCCCAGCGCTGCAGCCCCAGCCTGAGCAGCAGCCGTTCTCACAGATG
CTCCTGGGCCACCGTGGCCACATGCGTCGAGTGAAGACCACTACTCTGTAAATAGTACTGGCCGGCGG
GGGCGTGGCACCCTGGGACGGCTCCACCTGGACGGGGACGGAACCCAGGTGGGGGACCCCTGCGGCCT
GCAGCCTCCCTGCCTCACATTGCTAAGACTCAAAGGGATGCAGGCCATATTGCCAGCAAGAGCCCTGTC
ATGTTGGTGGCCCTGCGGCCAACCAACATGGACCGTGAAGGAGACAAGTTCTTCCAGTCCCATTAACAC
TACAATCCACAGTTTGAATACAGGAGCCCATGCCACGGCTGTGCTGGAAAAGTACTGCCAGGCCTCT
GGACAGTTCATCCATCAGGCAGTTGGCATCATTGAGGCTGTTCTGGAGAAGTTTGAACCTATGAACAC
TTTGAGGCTGCCACTGGGGGGCAGCTGCTACCAAGTGCAGATATGGTCGATTGTGCGCAAATACATG
CAGAAGGAGGGCTGCGTGGGGAGGTTGTGGTGCAGTGAAGGACCTGCTGTCCAGGCAGTGATG
ATGGTGGAGAACAGCCGGCCACATTGGCGATCAACCTGACCGGAGCCCGCCAGTACTGTTGGAGGGC
ATGCTGCGGCATGAGTAGGCACCACTACCTGCGGGGCGTGAACAACGCGCGCCAGCCGTGGCACAAC
GCGGAGGGCCGGCTGCGGTACGGGCTGCGGCCGCGGAACCCACGAGGAGGGCCTGGCCAGCCTGCAC
AGCGTGCTGTTCCGCAAGCAGCCGTTCCGTGGCGCGCTGCACTGCTCTACTACACCATCCACCGGCC
GCGCGCATGTCTTCCGTGAGCTCTTCCAGGACCTGGAGCGCTACGTGCAGGACGCCGATGTGCGCTGG
GAGTACTGCGTGCGGCCAAGCGCGGCCAGACCACTCGCTGCCAGGTTGTTTCAGCAAGGACCAG
GTGTACCTGGACGGCATCGTGCGATTCTGCGACATCGCCAGACCATCGATTTCCCGTTGCTGACCTCA
CTGGGCAAGGTGCTATGAGGATGTGGACCACTGCGGCCCCATGGGGTGTGGATAATACCCGGGTG
CCCCATTATGCAGGACTTGGCAGCTACCGGCAGCAGCTGGAGCACATCATGGCCACCAACCGGCTG
GATGAGGCGGAGCTGGTTCGCTGCTACCCGAC TGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCG



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001040715
Insert Size:	1416 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001040715.1</u>
RefSeq Size:	3416 bp
RefSeq ORF:	1416 bp
Locus ID:	653319
UniProt ID:	<u>Q68EN5</u>
Cytogenetics:	16q22.1
MW:	53.4 kDa