

Product datasheet for **SC311164**

UEVLD (NM_001040697) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	UEVLD
Synonyms:	ATTP; UEV3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF:

>SC311164 representing NM_001040697.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAAATTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGTTCGACTGCGAGGGCCTGAGACGGCTGCTTGGCAAGTACAAGTTCAGGGACCTAACTGTGGAA
GAACTAAGGAATGTAATGTATTTTCCACATTTCAATATTCATGGACACCTATGTTTTAAAGAT
AGTTCTCAGAAAGACCTGCTGAATTTACTGGCACAATTCCTGTGATGTATCAGGGTAATACATATAAC
ATACCAATTCGTTTCTGGATTTTGGATTCTCACCTTTTCGCTCCCCTATTGCTTCTTGAAGCCAACT
GCAAAATATGGGAATCTTAGTCGGAACATGTGGATGCTCAAGGCAGAATATTTGCCCTATCTCCAA
AACTGGAGCCATCTAAATCTGTCTATTGTTGGATTAATTAAGAAATGATTGCCAAGTTTCAAGAGGAA
CTTCCCATGTATTCTCTATCATCTGATGAGGCACGGCAGGTAGACTTGTAGCCTATATTGCAAAA
ATCACTGAAGGTGTTTCAGATACAAATCAAAGAGCTGGGCAATCATGAGAATAAAACAGTCAATAAA
ATTACTGTGGTTGGAGGTGGAGAACTCGGTATTGCCTGCACATTAGCAATTTCAAGAAAGGTATTGCA
GACAGGCTTGTCTCTTAGACCTCTCAGAAGGGACTAAAGGAGCCACGATGGACCTTGAAATCTTCAAC
CTTCTCAATGTGGAGATCAGCAAGATTTGTCTGCCTCTGCTCATTCCAAGGTGGTGTCTTCACAGTC
AACTCTTTGGGTAGTTCTCAGTCGTACCTTGATGTGGTACAGAGCAATGTGGATATGTTAGAGCCCTT
GTCCCAGCTCTGGGACATTATAGTCAACACAGTGTCTGCTCGTTGCATCTCAACCAAGTGGAAATCATG
ACCTATGTAAATGGAAGTGAATACATTTCTGCAATCGAGTATCGGAATTGGATGTAATCTGGAT
TCACAGAGATTACAGTATATTATTACAAATGTTTTGAAGGCACAGACTTCAGGCAAGAAAGTATGGGTT
ATTGGCGAGCAAGGAGAAGACAAAGTGTCTACATGGAGTGGCAAGAAGAAGTAGTGAGTCATACCTCT
CAAGTGCAGCTGTCCAACAGAGCCATGGAAGTCTAAGAGTAAAGGTCAAAGATCCTGGTCTGTTGGA
CTATCAGTAGCTGACATGGTTGACAGTATTGTAACAATAAGAAGAAAGTGCATTCTGTATCAGCTTTA
GCAAAAGGATATTATGATATAAATAGTGAAGTGTTTTAAGTTTGCCTTGCATCCTTGGAACCAATGGA
GTATCTGAAGTTATCAAACACACTGAAAGAAGATACAGTTACTGAGAACTCCAAGCAGTGCATCC
TCAATCCACAGTCTCCAACAACAGTTAAACTTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGCCCGC
```

Restriction Sites:

SgfI-MluI

ACCN:

NM_001040697

Insert Size:

1416 bp

OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:

This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001040697.2</u>
RefSeq Size:	2425 bp
RefSeq ORF:	1416 bp
Locus ID:	55293
UniProt ID:	<u>Q8IX04</u>
Cytogenetics:	11p15.1
MW:	52.3 kDa
Gene Summary:	Possible negative regulator of polyubiquitination.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longest isoform (α). Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.