

Product datasheet for **SC311221**

Galectin 7 (LGALS7B) (NM_001042507) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Galectin 7 (LGALS7B) (NM_001042507) Human Untagged Clone
Tag:	Tag Free
Symbol:	Galectin 7
Synonyms:	Gal-7; GAL7; HKL-14; LGALS7; PI7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC (PS100020)
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_001042507 edited ACGGCTGCCCAACCCGGTCCCAGCCATGTCCAACGTCCCCACAAGTCCCTCACTGCCCGA GGGCATCCGCCCTGGCACGGTGCTGAGAATTCGCGGCTTGTTCTCCCAATGCCAGCAG GTTCCATGTAAACCTGCTGTGCGGGGAGGAGCAGGGCTCCGATGCGGCCCTGCATTTCAA CCCCCGGCTGGACACGTGGAGGTGGTCTTCAACAGCAAGGAGCAAGGCTCCTGGGGCCG CGAGGAGCGCGGGCCGGCGTTCCTTTCCAGCGGGCAGCCCTTCGAGGTGCTCATCAT CGCGTCAGACGACGGCTTCAAGGCCGTGGTTGGGGACGCCAGTACCACCACTTCCGCCA CCGCTGCCGCTGGCGCGCTGCGCCTGGTGGAGGTGGGGGGGACGTGCAGCTGGACTC CGTGAGGATCTTCTGAGCAGAAGCCAGGCGGGCCCGGGGCTTGGCTGGCAAATAAAGC GTTAGCCCGCAGCGCAAAAAAAAAAAAAAAAAAAAA
Restriction Sites:	Please inquire
ACCN:	NM_001042507
Insert Size:	600 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	The ORF matches with reference, NM_001042507.1 except one SNP.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001042507.1](#), [NP_001035972.1](#)

RefSeq Size: 411 bp

RefSeq ORF: 411 bp

Locus ID: 653499

UniProt ID: [P47929](#)

Cytogenetics: 19q13.2

Gene Summary: The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. Differential and in situ hybridization studies indicate that this lectin is specifically expressed in keratinocytes and found mainly in stratified squamous epithelium. A duplicate copy of this gene (GeneID:3963) is found adjacent to, but on the opposite strand on chromosome 19. [provided by RefSeq, Jul 2008]