

## Product datasheet for SC311225

### CARHSP1 (NM\_001042476) Human Untagged Clone

#### Product data:


Product Type:	Expression Plasmids
Product Name:	CARHSP1 (NM_001042476) Human Untagged Clone
Tag:	Tag Free
Symbol:	CARHSP1
Synonyms:	CRHSP-24; CRHSP24; CSDC1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC311225 representing NM_001042476. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCATCTGAGCCTCCCCACCACCACAGCCCCCACCATCAAGTTTCAGTCGGGCTGCTGGACACC
CCTCGGAGCCGTGAGCGCTCACCATCCCCTCTGCCGGGCAACGTGGTCCCAAGCCCACTGCCCACTCGC
CGGACGAGGACCTTCTCGGCGACGGTGC GGCTTCACAGGGCCCCGTCTACAAAGGAGTCTGCAAATGC
TTCTGCCGGTCCAAGGGCCATGGCTTCATTACCCAGCTGATGGCGGCCCGACATCTTCTGCACATC
TCTGATGTGGAAGGGGAGTATGTCCAGTGAAGGCGACGAGTACCTATAAAATGTGCTCCATCCCA
CCCAAGAATGAGAAGCTGCAGGCCGTGGAGTCTGTCATCACTCACCTGGCACCAGGCCAAGCATGAG
ACCTGGTCTGGACATGTCATCAGCTCCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

**Plasmid Map:**


<b>ACCN:</b>	NM_001042476
<b>Insert Size:</b>	444 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001042476.2</a>
<b>RefSeq Size:</b>	3043 bp
<b>RefSeq ORF:</b>	444 bp
<b>Locus ID:</b>	23589
<b>UniProt ID:</b>	<a href="#">Q9Y2V2</a>
<b>Cytogenetics:</b>	16p13.2
<b>Protein Families:</b>	Transcription Factors
<b>MW:</b>	15.9 kDa
<b>Gene Summary:</b>	<p>Binds mRNA and regulates the stability of target mRNA. Binds single-stranded DNA (in vitro). [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1-9 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>