

Product datasheet for **SC311254**

FAM98B (NM_001042429) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM98B (NM_001042429) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAM98B
Synonyms:	FLJ38426
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC311254 representing NM_001042429. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
ATGAGAGGGCCGGAGCCGGGTCCCAACCGACGATGGAGGGAGACGTGCTGGACACACTGGAGGCGCTG
GGGTATAAAGGACCATTGTTAGAAGAGCAAGCCCTTACAAAGCGGCAGAGGGTGGATTATCTTACCT
GAATTTTCAGAGCTCTGTATTTGGTTAGGCTCTCAAATAAAATCATTATGCAACTTGAAGAAAGTATC
ACGCTGCTGGAAGAGATGATCTAGAGAGCTTCCAGCTTGAGATAAGTGGCTTTTTAAAGAAATGGCA
TGTCATATTCTGTACTCATATCAGGAGATATTAAGATCGTTTAAAAAGAAGGAGGACTGTTGAAA
CTTCTATTATTTTAAGTACAGAACTTCAAGCTTCACAGATATTACAGAACAAAGAAACATAAAATTTCT
CAATTAGATAAAATAGTGAAGTTTATCAGGAAGTTCAAGCTATGTTTGATACACTTGGTATACCCAAG
TCAACAACCTTCTGACATTCGCGATATGCTAAACCAAGTGAATCAAAGGTGAAAGATATTCTCTCAAAG
GTCCAGAAAAATCATGTGGGAAAACCACTGCTGAAAATGGATTTAAATTCAGAACAGGCGGAACCACTG
GAAAGAATCAATGATGCTCTTTCCTGTGAATATGAGTGCCGCCGACGAATGTTAATGAAACGATTAGAT
GTGACTGTACAGTCCTTTGGATGGTCTGATAGAGCAAAGTAAAAACAGATGATATAGCAAGAATTTAT
CAGCCTAAGCGTTATGCTTTGTACCCAAGACAACGATTACAATGGCACATCTACTTGCTGCTCGTGAA
GATCTATCCAAGATCATTAGGACAAGTAGTGGCACCAGCCGGGAGAAGACCGCATGTGCCATTAATAAG
GTTGGTGTTCCTTCAGTACAGTGAAAAATGAATTGATGATATCTTACTTAATGTTTTTACAGATTTTA
GTCTATTTTCTTTCATGTCTGGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001042429
Insert Size:	993 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001042429.1</u>
RefSeq Size:	1101 bp
RefSeq ORF:	993 bp
Locus ID:	283742
Cytogenetics:	15q14
MW:	37.2 kDa
Gene Summary:	Positively stimulates PRMT1-induced protein arginine dimethylated arginine methylation (PubMed:28040436). Promotes colorectal cancer cell malignancy (PubMed:28040436). [UniProtKB/Swiss-Prot Function]