

## Product datasheet for **SC311281**

### ACBD5 (NM\_001042473) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ACBD5 (NM_001042473) Human Untagged Clone
Tag:	Tag Free
Symbol:	ACBD5
Synonyms:	RDLKD
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC311281 representing NM\_001042473.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCGGACACGAGATCCGTGCACGAGACTAGGTTTGAGGCGCCGTGAAGGTGATCCAGAGTTTGCCG
AAGAATGGTTTCATTCCAGCCAACAAATGAAATGATGCTTAAATTTATAGCTTCTATAAGCAGGCAACT
GAAGGACCTGTAAACTTTCAAGGCCTGGATTTTGGGATCCTATTGGAAGATATAAATGGGATGCTTGG
AGTTCACTGGGTGATATGACCAAAGAGGAAGCCATGATTGCATATGTTGAAGAAATGAAAAAGATTATT
GAAACTATGCCAATGACTGAGAAAGTTGAAGAATTGCTGCGTGCATAGGTCCATTTTATGAAATTGTC
GAGGACAAAAAGAGTGGCAGGAGTTCTGATATAACCTCAGATCTTGGTAATGTTCTCACTTCTACTCCA
AACGCCAAAACCGTTAATGGTAAAGCTGAAAGCAGTGACAGTGGAGCCGAGTCTGAGGAAGAAGAGGCC
CAAGAAGAAGTAAAGGAGCAGAACAAGTGATAATGATAAGAAAATGATGAAGAAGTCAGCAGACCAT
AAGAATTTGGAAGTCATTGTCACTAATGGCTATGATAAAGATGGCTTTGTTCCAGGATATACAGAATGAC
ATTCATGCCAGTCTTCCCTGAATGGCAGAAGCACTGAAGAAGTAAAGCCATTGATGAAAACCTGGGG
CAAACCTGGAAAATCTGCTGTTTGCATTCACCAAGATATAAATGATGATCATGTTGAAGATGTTACAGGA
ATTCAGCATTTGACAAGCGATTGACAGTGAAGTTTACTGTGATTCTATGGAACAATTTGGACAAGAA
GAGTCTTTAGACAGCTTTACGTCCAACAATGGACCATTTCACTTACTTGGGTGGTCAATCCAGTCAA
CCCATGGAAAATCTGGATTTCTGTGAAGATATTCAAGTACCTCCTGGAAATGGCAACATTGGGAATATG
CAGGTGGTTGCAGTTGAAGGAAAAGGTGAAGTCAAGCATGGAGGAGAAGATGGCAGGAATAACAGCGGA
GCACCACACCGGGAGAAGCGAGGCGGAGAACTGACGAATCTCTAATGTTAGAAGAGGAAGAGGCAT
AGGATGCAACACTTGAGCGAAGGAACCAAGGGCCGGCAGGTGGGAAGTGGAGGTGATGGGAGCCGCTGG
GGCTCCGACAGAGGGTCCCGAGGCAGCCTCAATGAGCAGATCGCCCTCGTGTGATGAGACTGCAGGAG
GACATGCAGAATGTCCTTCAGAGACTGCAGAAACTGGAACGCTGACTGCTTTGCAGGCAAAATCATCA
ACATCAACATTGCAGACTGCTCCTCAGCCACCTCACAGAGACCATCTTGGTGGCCCTTCGAGATGTCT
CCTGGTGTGCTAACGTTTGCCATCATATGGCCTTTATTGCACAGTGGTGGTATTATATACTATCAA
AGAAGGAGAAGAAAACTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:**

SgfI-MluI

**Plasmid Map:**



**ACCN:**

NM\_001042473

<b>Insert Size:</b>	1473 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001042473.3</a>
<b>RefSeq Size:</b>	3871 bp
<b>RefSeq ORF:</b>	1473 bp
<b>Locus ID:</b>	91452
<b>UniProt ID:</b>	<a href="#">Q5T8D3</a>
<b>Cytogenetics:</b>	10p12.1
<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	54.8 kDa
<b>Gene Summary:</b>	<p>This gene encodes a member of the acyl-Coenzyme A binding protein family, known to function in the transport and distribution of long chain acyl-Coenzyme A in cells. This gene may play a role in the differentiation of megakaryocytes and formation of platelets. A related protein in yeast is involved in autophagy of peroxisomes. A mutation in this gene has been associated with autosomal dominant thrombocytopenia. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR, lacks part of the 5' coding region, and uses a downstream start codon, compared to variant 1. The encoded isoform (2) has a shorter N-terminus, compared to isoform 1. Variants 2 and 17-19 all encode the same isoform (2).</p>