

Product datasheet for **SC311307**

ANKZF1 (NM_001042410) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ANKZF1 (NM_001042410) Human Untagged Clone
Tag:	Tag Free
Symbol:	ANKZF1
Synonyms:	Vms1; ZNF744
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC311307 representing NM_001042410.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCGCCGGCTCCAGATGCAGCCCCGGCTCCTGCGTCGATCTCCCTGTTTGACCTCAGCGCGGATGCT
CCGGTCTTTTCAGGGCTGAGCCTGGTGAGCCACGCGCTGGGAGGCTCTGGCCCCGGCTCCGCGTACT
TCCTGTTTCAGGCTCAGGGGAGAGAGAAAGCCAGAAAGAAAGCTACTCCAGGGTCTATGGATATTTCA
GAGAAGTTATTTTGTCAACTTGTGACACAGCCTTCCAGAACCACCAAGAAGACAGAGGGAACATTATAAG
CTTGACTGGCATCGTTTAACTAAAGCAACGTCTCAAGGACAAGCCTCTCCTGTCTGCCCTGGACTTT
GAAAAGCAGAGCTCCACAGGAGATCTTTCCAGCATCTCGGGATCAGAAGACTCAGACTCAGCCAGTGAG
GAGGACTTGCAGACACTGGATCGGGAGAGGGCTACATTTGAGAAGTTGAGCCGACCCCCAGGCTTTTAC
CCTCATCGAGTTCTTTCCAGAATGCCAGGGCCAGTTTCTTTATGCCTACCGCTGTGCTCAGGCCCT
CATCAGGATCCCCAGAAGAGGCAGAAGTGTCTACAGAACCTGCAAAGTAGAGGTCCCAGAGACTGC
GTGGTGCTCATGGCTGCAGCTGGGCACTTGTGTTGCTATATTTCAAGGAAGAGAAGTGGTGACACAC
AAAATTTTACCCTATACGGTTCGGGCAAGCGGGGACAGCCAGGGGCTTCGGGATGCCCGAGGT
GGGCCATCACACTCTGCTGGAGCAACCTGAGGCGCTACAATGAAGCCACACTATAAAGGATGTTGCT
GACCTGCTGGCAGGGCAAGCTGGGTAAGGCGCTGGAGGAGGCTGGTACAATACTGTTGCGTGCTCCC
CGCTCTGGCCGGTCTTTGTTCTTTGGAGGCAAGGGAGCACCCCTGCAAAGGGGGGATCCCCGACTTTGG
GATATCCCCCTCGCCACCCGAGACCCACCTTCCAAGAGCTACAGCGTGTGCTCCATAAGCTGACCACT
TTGCATGTCTATGAAGAAGACCTCGGGAAGCAGTCAGACTGCACTCACCTCAGACACTGGAAAAACA
GTAAAGAGAGGAGAGAAAAGCCTACTGAGGAAGAAAATAAGAAAGATCTGCAGGGATGAAAAGGAAGCG
CTGGGGCAGAATGAGGAATCTCCAAACAGGGTTCAGGGTCCGAGGGGAGAAGATGGCTTTTCAGGTAGAG
TTGGAGCTAGTGGAGTTGACTGTGGGACTCTGGATCTTTGTGAGTCTGAAGTATTGCCCAAGCGGAGG
AGGAGAAAAAGGAATAAGAAGGAGAAAAGCCGAGACCAGGAGGCTGGGGCACATCGGACTCTTCTCCAG
CAAACAAGAAGAGGAGCCTTCCACACAGTCATCCAGGCAGTTGCTGCCCCCTTGGGCCCTTTGCTG
GATGAGGCCAAAGCCCCTGGTCAGCCAGAGCTCTGGAATGCACTGCTTGTGCTTGGCAGCTGGAGAT
GTTGGAGTGCTAAAGCTGCAGCTAGCTCCAGCCCTGCAGACCCTAGAGTTCTGTCTCTGCTCAGTGCC
CCCTTGGGCTCCGGTGGCTTTACTCTCTGCATGCAGCAGCTGCAGCTGGAAGAGGCTCAGTGGTTCGT
CTGCTGCTGGAAGCAGGTGCTGACCCACTGTGCAGGACTCTCGGGCCCCGGCCACCTTATACTGTTGGC
GCTGACAAATCAACACGTAAAGTTCGAAGGTTTCATGGAGAAGAAATCCAGATGCCTACGATTACAAC
AAGGCTCAGGTGCCAGGACCATTTGACACCAGAAAATGGAGGCACGGCAGGCTACACGGAAAAGGGAGCAG
AAGGCAGCCCCGGCGCAACGGGAGGAACAGCAGAGGAGCAGGAGCAGGAGGAGCGTGAACGAGAA
GAGCAGCGGCGATTTGCCGCCCTCAGTGACCGAGAGAAGAGAGCTCTGGCTGCAGAGCGCCGACTCGCT
GCCAGTTGGGAGCCCCTACCTCTCCAATCCCTGACTCTGCAATCGTCAATACTCGACGCTGCTGGAGT
TGTGGGGCATCCCTCAAAGGCCTGACTCCCTTTCACTACCTCGACTTCTCTTTCTGCTCCACACGTTGC
CTCCAGGATCATCGCCGTCAAGCAGGGAGGCCCTTTCCGTA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001042410

Insert Size: 2181 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001042410.1](#)

RefSeq Size: 2521 bp

RefSeq ORF: 2181 bp

Locus ID: 55139

UniProt ID: [Q9H8Y5](#)

Cytogenetics: 2q35

MW: 80.9 kDa

Gene Summary: Plays a role in the cellular response to hydrogen peroxide and in the maintenance of mitochondrial integrity under conditions of cellular stress (PubMed:28302725). Involved in the endoplasmic reticulum (ER)-associated degradation (ERAD) pathway (By similarity).
[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (2) has an alternate 5' UTR exon, compared to variant 1. Variants 1 and 2 encode the same isoform (1).