

## Product datasheet for **SC311311**

### Calpastatin (CAST) (NM\_001042441) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Calpastatin (CAST) (NM_001042441) Human Untagged Clone
Tag:	Tag Free
Symbol:	CAST
Synonyms:	BS-17; PLACK
Vector:	<u><a href="#">pCMV6 series</a></u>



[View online »](#)

**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_001042441, the custom clone sequence may differ by one or more nucleotides

```

ATGTCCCAGCCCGGCCAGAAAGCCCGCCGCTCCCCCGGGCCCCGGCGAGCAGCCGCCGCC
CGCCGCACCCATGAGCATGTCAGTGAAAAAACAGTGAATCGCCTTCCAAACCAGGAGAA
AAGAAAGGATCAGATGAGAAAAAGCAGCAAGCCTCGGCAGCAGTCAATCCTCCAGAACC
TATGCTGGTGGAACAGCCTCGGCCACCAAGGTGTCAGCTTCTCTGGTGAACCAGCAAG
TCTTCCAGTATGAATCCACAGAAACCAAGGCCATTCCAGTCAGCCAACAGATGGAAGGA
CCACATCTTCTAACAAAGAAAAACAAAAAACAGGCTGTAAAAACAGAACCTGAGAAG
AAGTCACAGTCAACCAAGCCAAAAAGCCTACCCAAGCAGGCATCAGATACAGGAAGTAAC
GATGCTCACAAATAAAAAAGCAGTTTCCAGATCAGCTGAACAGCAGCCATCAGAGAAATCA
ACAGAACCAAAGACTAAACCACAAGACATGATTTCTGCTGGTGGAGAGAGTGTGCTGGT
ATCACTGCAATATCTGGCAAGCCGGGTGACAAGAAAAAGAAAAGAAATCATTAAACCCCA
GCTGTGCCAGTTGAATCTAAACCGGATAAACCATCGGGAAAGTCAGGCATGGATGCTGCT
TTGGATGACTTAATAGATACTTTAGGAGGACCTGAAGAACTGAAGAAGAAAATACAACG
TATACTGGACCAGAAGTTTCAGATCCAATGAGTTCACCTACATAGAGGAATTGGGTAAA
AGAGAAGTCACAATTCTCCAAAATATAGGGAACATTGGCTAAAAAGGAAGGGATCACA
GGGCTCTGCAGACTCTTCGAAACCCATAGGGCCAGATGATGCTATAGACGCCTTGTCAT
TCTGACTTACCTGTGGGTGCCTACAGCTGCTGGAAAGAAAACTGAAAAAGAGGAATCT
ACAGAAGTTTTAAAGCTCAGTCAGCAGGGACAGTCAGAAGTGTGCTCCACCCCAAGAG
AAGAAAAGAAAGGTGGAGAAGGATACAATGAGTGTCAAGCACTCGAGGCTCTGTGGCT
TCACTGGGCACCCGGCAAGCAGAACCTGAGCTCGACCTCCGCTCAATTAAGGAAGTCGAT
GAGGCAAAAGCTAAAGAAGAAAACTAGAGAAGTGTGGTGAAGGATGATGAAACAATCCCA
TCTGAGTACAGATTAACCAAGCCACGGATAAAGATGGAACCACTATTGCCAGAGCCT
GAAGAAAAACCCAAAGCCTCGGAGTGAATCAGAACTCATTGATGAACCTTTCAGAAGATTTT
GACCGGTCTGAATGTAAGAGAAACCATCTAAGCCAAGTAAAAAGACAGAAGAATCTAAG
GCCGCTGCTCCAGCTCTGTGTCGGAGGCTGTGTGTCGGACCTCCATGTGTAGTATACAG
TCAGCACCCCTGAGCCGGCTACCTTGAAGGGCACAGTCCAGATGATGCTGTAGAAGCC
TTGGCTGATAGCCTGGGAAAAAGGAAGCAGATCCAGAAGATGGAACCTGTGATGGAT
AAAGTCAAGGAGAAGGCCAAAGAAGAAGACCGTAAAAAGCTTGGTAAAAAGAAGAAACA
ATTCTCTGATTATAGATTAGAAGAGGTCAAGGATAAAGATGGAAGGCCACTCTGCCA
AAAGAGTCTAAGGAACAGCTTCCACCCATGAGTGAAGACTTCTTCTGGATGCTTTGTCT
GAGGACTTCTCTGGTCCAAAAATGCTTATCTCTTAAATTTGAAGATGCTAAACTTGCT
GCTGCCATCTCTGAAGTGGTTTTCCAAACCCAGCTTCAACGACCCAAGCTGGAGCCCCA
CCCCGTGATACCTCGCAGAGTGACAAAGACCTCGATGATGCCTTGGATAAACTCTCTGAC
AGTCTAGGACAAAAGGCAGCCTGACCCAGATGAGAACAAACCAATGGAAGATAAAGTAAAG
GAAAAAGCTAAAGCTGAACATAGAGACAAGCTTGGAGAAAAGAGATGACACTATCCACCT
GAATACAGACATCTCTGGATGATAATGGACAGGACAAACCAAGTGAAGCCACCTACAAG
AAATCAGAGGATCAAGAAAACCTGCAGATGACCAAGACCCCATGATGCTCTCTCAGGA
GATCTGGACAGCTGTCCCTCCACTACAGAAACCTCACAGAACACAGCAAAGGATAAAGTGC
AAGAAGGCTGCTTCCAGCTCCAAAGCACCTAAGAATGGAGGTAAGCGAAGGATTCAGCA
AAGACAACAGAGGAACTTCCAAGCCAAAAGATGACTAA

```

**Restriction Sites:** Please inquire

**ACCN:** NM\_001042441

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001042441.1</a> , <a href="#">NP_001035906.1</a>
<b>RefSeq Size:</b>	4572 bp
<b>RefSeq ORF:</b>	2319 bp
<b>Locus ID:</b>	831
<b>UniProt ID:</b>	<a href="#">P20810</a>
<b>Cytogenetics:</b>	5q15
<b>Gene Summary:</b>	<p>The protein encoded by this gene is an endogenous calpain (calcium-dependent cysteine protease) inhibitor. It consists of an N-terminal domain L and four repetitive calpain-inhibition domains (domains 1-4), and it is involved in the proteolysis of amyloid precursor protein. The calpain/calpastatin system is involved in numerous membrane fusion events, such as neural vesicle exocytosis and platelet and red-cell aggregation. The encoded protein is also thought to affect the expression levels of genes encoding structural or regulatory proteins. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jun 2010]</p>