

Product datasheet for **SC311333**

ZNF706 (NM_001042510) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF706 (NM_001042510) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF706
Synonyms:	HSPC038; PNAS-106; PNAS-113
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC311333 representing NM_001042510. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
ATGGCTCGTGGACAGCAGAAAATTCAGTCTCAGCAGAAAATGCCAAAAGCAAGCTGGACAAAAGAAG
AAACAAGGACATGACCAAAAGGCTGCTGCCAAAGCTGCCTTAATATACCTGCACTGTCTGTAGGACA
CAATGCCAGACCCTAAGACCTCAAGCAGCACTTTGAGAGCAAGCATCCTAAGACTCCACTTCTCCA
GAATTAGCTGATGTTCAAGCAT**AA**
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001042510
Insert Size:	231 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001042510.1</u>
RefSeq Size:	2870 bp
RefSeq ORF:	231 bp
Locus ID:	51123
UniProt ID:	<u>Q9Y5V0</u>
Cytogenetics:	8q22.3
MW:	8.5 kDa
Gene Summary:	<p>Transcription repressor involved in the exit of embryonic stem cells (ESCs) from self-renewal. Acts by repressing expression of KLF4.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) encodes the functional protein. Variants 1, 3, 4 and 5 encode the same protein.</p>