

## Product datasheet for SC311334

### ZNF706 (NM\_001042511) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF706 (NM_001042511) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF706
Synonyms:	HSPC038; PNAS-106; PNAS-113
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC311334 representing NM_001042511. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG  
 GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**  
 ATGGCTCGTGGACAGCAGAAAATTCAGTCTCAGCAGAAAAATGCCAAAAGCAAGCTGGACAAAAGAAG  
 AAACAAGGACATGACCAAAAGGCTGCTGCCAAAGCTGCCTTAATATACCTGCACTGTCTGTAGGACA  
 CAAATGCCAGACCCTAAGACCTCAAGCAGCACTTTGAGAGCAAGCATCCTAAGACTCCACTTCTCCA  
 GAATTAGCTGATGTTCAAGCAT**AA**  
**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT  
 TACAAGGATGACGACGATAAGGTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001042511
Insert Size:	231 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.


[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001042511.1</a></u>
<b>RefSeq Size:</b>	2673 bp
<b>RefSeq ORF:</b>	231 bp
<b>Locus ID:</b>	51123
<b>Cytogenetics:</b>	8q22.3
<b>MW:</b>	8.5 kDa
<b>Gene Summary:</b>	Transcription repressor involved in the exit of embryonic stem cells (ESCs) from self-renewal. Acts by repressing expression of KLF4.[UniProtKB/Swiss-Prot Function]