

Product datasheet for SC311379

TTC14 (NM_001042601) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TTC14 (NM_001042601) Human Untagged Clone
Tag:	Tag Free
Symbol:	TTC14
Synonyms:	DRDL5813; PRO19630
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC311379 representing NM_001042601. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGACCGGGACCTTTTGCAGCAGTCGCTAAATTGCCACGGGTCGCTTTGCTCTCTCTACTTCGGAGC
GAACAGCAGGACAATCCACACTTCCGTAGCCTCCTGGGGTCGGCCGCCGAGCCAGCCGGGGCCCGCCG
CCCCAGCACCCGTTGCAGGGCAGAAAAGAGAAGAGAGTTGACAACATCGAGATACAGAAATTCATCTCC
AAAAAAGCGGATCTGCTTTTGCACCTTCTGGAAATCAGATGCACCTGCAACTTCTGAAATTAATGAA
GACAGTGAAGATCATTATGCAATCATGCCACCTTATAGCAATTCATGGAGATACCTAGTATGGATCGG
AGAGAGCTGTTTTCCGAGATATTGAGCGTGGTGATATAGTGATTGGAAGAATTAGTTCTATTCCGGAA
TTCGGTTTTTTCATGGTGTGATCTGTTAGGAAGTGGTATCATGAGAGATATAGCCCACTAGAAATC
ACAGCTCTTTGTCCCTTAAGAGATGTGCCCTTCTCACAGTAACCATGGGGATCCTTTATCATATTACCAA
ACTGGTGACATCATTCGAGCTGGAATCAAGGATATTGACAGATACCATGAAAAGCTAGCAGTATCTCTG
TATAGCTCTTCTCTCCACCACACCTATCTGGTATTAAATTAGGTGTAATTAGCTCTGAAGAGCTTCTT
TTATACTACAGGAGAAGTGTGAGCTAAATAGCAATTCCTTGGAGTCCTATGAAAATGTCATGCAGAGT
TCCTTGGGATTTGTTAATCCAGGAGTAGTTGAATTCCTTCTAGAAAACTAGGAATAGATGAATCTAAT
CCACCATCTTTAATGAGAGGCCACAAAAGCAAAAATTTCTCTGAAGATGATTTTGCTTCTGCATTGAGA
AAAAACAATCCGCATCTTGGGCTTTAAAATGTGTGAAGATCGGAGTTGACTATTTTAAAGTTGGACGC
CATGTGGATGCTATGAATGAATACAATAAAGCTTTGGAATAGACAAAACAAACGTTGGAAGCTTTGGTA
GCTCGTGGAGCATTATATGCGACAAAAGGAAGTTTGAACAAAGCAATAGAAGATTTTGAGCTTGCAATTA
GAAAATGTCCAACACACAGAAATGCAAGAAAATACCTCTGCCAGACACTGTAGAGAGAGGAGGACAG
TTAGAAGAAGAAGAAAAGTTTTTAAATGCTGAAAGTTACTATAAGAAAGCTTTGGCTTTGGATGAGACT
TTTAAAGATGCAGAGGATGCTTTCAGAAACTTCATAAATATATGCAGGTGATTCCTTATTTCTCTTA
GAAATTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001042601
Insert Size:	1320 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001042601.2</u>
RefSeq Size:	4813 bp
RefSeq ORF:	1320 bp
Locus ID:	151613
UniProt ID:	<u>Q96N46</u>
Cytogenetics:	3q26.33
MW:	49.7 kDa