

Product datasheet for **SC311402**

ZNF699 (NM_198535) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF699 (NM_198535) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF699
Synonyms:	hang
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC311402 representing NM_198535.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGGAAGAAAGAAAACTGCTGAGTTACAGAAAAATAGAATACAGGACTCAGTAGTCTTTGAGGAT
GTGGCTGTGGACTTTACCCAGGAGGAATGGGCTTTGCTGGATCTTGCTCAGAGAAACCTCTACAGAGAT
GTGATGCTGGAAAACTTCCAGAACTGGCCTCACTAGGTATCCGCTACACACACCCCATCTGATCTCC
CAGTGGGAACAAGAGGAGGACCTGCAGACAGTGAAGAGAGAACTTATCCAGGGCATCTTTATGGGAGAG
CACCAGGAAGGTTTTGAGACTCAACTTAAACTAATGAATCAGTTGCTTCACAAGATATTTGTGGAGAA
AAAATATCCAATGAACAGAAAAATAGTAAGATTCAAAAGGAATGATTCTGGTTCTCCAGTTTACATGAA
AACCAGGAGTCTGTGGTATTGATTATCAGAACAAAAGCCATGAGAGACATTTGAGAAATCATATGGTA
GAGAATCTATGAATGTTATGAAGAAAATCAAGATGGACAGACTTTAGCCAAGTTCCAAATCTTGAT
TCACTCAAGAGAAATACTGAAGTAAATCCTGTGAATGCCATGAGTGTGGAAGGCCCTTCGTGGATCAT
TCATCCCTTAAGAGTCACATCAGGTCTCATACTGGAAGCAAACCTATCAGTGCAAGGAATGTGGGAAG
GCCTTCCATTTTCTGCTTGTTCAGAAAGCATATGAAAACCCCACTGAAGAGAAGCCCTATGAATGT
AAGGAATGTACCAAGCCTTCAGCTGTTCTCATTCTTTAGGGCAGATATGAAGATTCACATCGGAAG
ACAACTATGAATGTAAGAATGTGGGAAAGGTTTTAGTTGTTCTCATCGCTCAGAGACACAAAAAGA
ATTCACAGTGGAGATAAGCCTTATGAATGTAAGGAATGTGGGAAGGCCTTCAGTTGTTCTCCTCACTC
TCCAAACACAAAAGAAATTCACAGTGGAGATAAGCCATATGAATGTAAGGAATGTGGGAAGGCCTTAGC
TCTTCTCTCACCTTATAATACATATAAGAATTCACACTGGAGAGAAGCCTTATGAATGTAAGAGTGT
GGGAAGGCCTTCAGCGAGTCCTCAAACTCACTGTACATGGGAGGACACATACTGGAGAGAAACCTAT
AAATGAAGGAATGTGGGAAAGCCTACAATTGTCCTCCTCTTAAGTATCCATATGAGAAAACACAT
GGAGAAAAACCGTATGAATGTCTGGAATGTGGAAAGGCCTTCTACCTTCCCACTTCCCTTAATACACAT
GTGAAAAATCAAAGTCGAGAGAAACCTATGAATGTAAGGAATGTGGGAAGGCCTTTAGTTGTCCCTCG
TCCTTTAGAGCACATGTGAGAGATCACACTGGAAGATACAGTATGAATGTAAGGAGTGTGGGAAGACC
TTAGTCGTTCTCATCCCTCACCAGACACCTAAGAACTCACAGCGGAGAGAAGCCGTATGAATGTAAA
GAATGTGGGAAGCCTTTATTAGTTCTCCACCTTACAGTACATATCAGAACTCACACTGGAGAGAAA
CCCTATGAATGTAAGAAATGTGGGAAAGCCTTTATTTATCCCTCAGCCCTTAGGATCCACATGAGAACG
CACACTGGGAAAAACCTATGAATGTAAGGAATGTGGGAAAGCATTTTCGTATTCTTCATACCTTACT
GTACATGCAAGAAATGCACACTGGAGAGAAACCTTTGAATGTCTGGAATGTGGAAAGCTTTCAGTTGT
CCCTCATCTTTTGAAGGCATGTGAGAAGCCACACTGGAGAGAAACCTATGAATGTAAGGAATGTGGG
AAAGCCTTTGTTTGCCTGCTACTTTTGAAGGCATGTGAAAACCTCACACTAGAGAAAAACATA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_198535

Insert Size: 1929 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_198535.2</u>
RefSeq Size:	6634 bp
RefSeq ORF:	1929 bp
Locus ID:	374879
UniProt ID:	<u>Q32M78</u>
Cytogenetics:	19p13.2
MW:	74 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]