

## Product datasheet for **SC311407**

### ENGASE (NM\_001042573) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ENGASE (NM_001042573) Human Untagged Clone
Tag:	Tag Free
Symbol:	ENGASE
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC311407 representing NM\_001042573.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGAGCCCGCGCGGTGACGGTCACCCGTCGGCTACACGGCGCGCGCGGCAGCTGCAGGGGCTG
GCGGCCCGGAGCGGGGACGCAGGAGGAGCAGGAGGATCAGGAGCCCGGCCCGCGCGCGCGCGCG
GGAAGGAGCATCAAAGATGAAGAAGAAGAGACAGTCTTTTCGAGAGGTGGTCAATTTTTCCCGGACCCC
CTGCCAGTTAGATATTATGACAAGGACACCACAAACCAATCAGCTTTTACTTGTCTTCGTGGAGGAG
CTTTGGCGTGGAAAGCCCGCTTGGAGGATGGCTTAAATGTGGCCCTGGAGCCCTGGCGTGTGCCAG
CCCCCTGAGCAGCCAGAGGCCCGGACTTTGTTGTGCATGACATGATGGGCGGTACCTGGATGAC
AGGTTCAATCAGGGCTCGGTGGTGCAGACTCCCTATGCTTTCTACCACTGGCAGTGCATCGACGCTTTT
GTGTAATTCAGCCACCACACCGTACCATTCCCCAGTGGGCTGGACCAACTGCCACAGGCATGGG
GTCTGCGTGTGGGACTTTCATCACGGAGTGAATGAAGGGGAAGGCTCTGTGAAGCCTTCTGGCC
GGGATGAGCGCTCGTACCAGGCAGTGGCTGACCGGCTGGTCCAGATCACTCAGTTTTTTGTTTTGAT
GGCTGGCTGATCAACATCGAAGACTCGCTGAGTCTGGCCGCTGTGGGGAACATGCCTCCTTCTGCCG
TACCTCACCACACAGCTGCACCGCAGGTCCCAGGGGGCCTGGTCTCTGGTATGACAGCGTGGTCAA
AGTGGGAGCTCAAATGGCAAGACGAACTCAACCAGCACAAACAGGGTCTTCTTTGATTCTGCGACGGC
TTCTTCACTAACTATAACTGGCGGGAGGAGCACTTGGAGCGGATGCTGGGCGAGGCTGGGAGCGCCGG
GCTGATGTGTACGTGGGCGTGGATGTGTTTGTCTGAGGGAACGTGGTGGAGGCCGATTCGACACAGAC
AAGTCGTTGGAGCTGATCCGAAAGCATGGCTTCTCCGTGGCTTTGTTTGGCCCGGCTGGGTGTATGAG
TGCTGGAGAAGAAGGATTTCTCCAGAACCAGGACAAGTTCTGGGCGGACTGGAGCGTTATCTGCC
ACACATGACATCTGCTCCTTGCCTTTCTGTCAGTCTCTTCTGCTGGCATGGGTGCACGGAGGCTGTC
TATGGCCAGGAAGGCGGTAGGGCCCTGGTACCACCTGAGCGCCAGGAGATCCAGCCCTTGTGGGA
GAACACAGGCTGGGAGGGGATGGCCGGGCTGGGTGAGGACGCACTGCTGCCTGGAGGATGCCTGGCAC
GGAGGCAGCTCCCTGCTCGTCCGGGTGTGATCCACCGGAGGTTGAAATGTGGCTGTGAGGTTATTT
TCCCTGCAGGCCCCAGTGCCACCCAAGATTTACCTGTCCATGGTGTATAAGCTTGAGGGGCCACGGAC
GTCACAGTTGCTTTGGAGCTGACCACAGGGGATGCCGGCAGCTGCCACATCGGTGGCATCTCAGTGTG
AACGCAGAAACAAGCTCAAGACACAGCCTCCGACCCCTCCGGTGGCCCCACCAAGCTGGCCAGATGG
GTGGGCCGCTGCGGCCGGCAGCTGAGTGGGGGCTGGGTCCAGCACTGCTACGAGGTGAGCCTGCGTGGG
TGCTGCTGCTAGACCTCCTCGTTTGTCTCTACGGCCCGGGTGTAGTGGGAGGAGGAGCTTACC
TGTGCGCTTGGAGAGATCCAGGTGGTGGAGCTGCCAGCCTGCTGGCCCTCTGCCCCAGGTGCAGGCC
GTCACCATCTCTCACATCCGCTGGCAGCCATCCGCTCTGAGCGGGAGGGGCCCTGCTGCTCCAG
CTCAGTGCACCTGCACTGGTCTTCTCTCTCACAAGTCCGTTGCTTCCGAATCCACTGCTGGGGA
GGGATGAGTGTGACTCTCCGGGAGGGAGCTGCCAGGCCAGAGATGCCATGTTCTGGGTTGGCT
TTTGCCACCCAGTACCGGATAGTGGACCTGCTGGTGGAAAGCCCGGGCCCGGCCAGGATCGTCGCATG
GAATTTCTGGTGGAGCCTGTCCCAAGGAAGGTTCCGGGTACCTCAGGCCGAGTGGGCGAGGCAGTT
CTGCTTTATTCAGCCCTGCATGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** Sgfl-Mlul

**Plasmid Map:**


**ACCN:** NM\_001042573

**Insert Size:** 2232 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001042573.2](#)

**RefSeq Size:** 4494 bp

**RefSeq ORF:** 2232 bp

**Locus ID:** 64772

**UniProt ID:** [Q8NFI3](#)

**Cytogenetics:** 17q25.3

**Protein Pathways:** Other glycan degradation

**MW:** 84 kDa

**Gene Summary:** This gene encodes a cytosolic enzyme which catalyzes the hydrolysis of peptides and proteins with mannose modifications to produce free oligosaccharides. [provided by RefSeq, Feb 2012]