

## Product datasheet for **SC311409**

### ZNF33A (NM\_006954) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF33A (NM_006954) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF33A
Synonyms:	KOX2; KOX5; KOX31; NF11A; ZNF11; ZNF11A; ZNF33; ZZAPK
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


**Fully Sequenced ORF:** >SC311409 representing NM\_006954.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAACAAGGTAGAACAAGAAGTCCCAGGAGTCAGTATCATTAAAGATGTGACTGTGGGCTTCCCCAG
GAGGAGTGGCAGCACCTGGACCCTAGTCAGAGGGCTCTGTATAGAGATGTGATGTGGGAACTACAGC
AACCTTGTCTCAGTGGGTATTGTGTTCAAACCCAGAGGTGATCTTCAGGCTGCAACAAGGAGAAGAG
CCATGGAAACAGGAGGAAGAATTCCAAGCCAAAGCTTTCCAGAAGTCTGGACAGCTGATCACCTGAAA
GAGAGGAGCCAAGAAAACCAATCTAAACATTTGTGGGAAGTTGTATTCAATAAATGAAATGCTGACT
AAGGAACAAGGTGATGTAATAGGAATACCATTTAATGTGGATGTAAGTTCTTTCCCTCCAGAAAAATG
TTCTGTCAAGTGTGATTGATGTAATAGGATGAGTTTCAACACTGTTTCCAGAATTGGTTATCAGTAAAGATAAAC
TATTTAGGAAAAAAGTCTGATGAATTTAATGCCTGTGGGAAATGTTACTCAATATTAAGCATGATGAA
ACTCATACTCAAGAGAAAAATGAAGTTTGAAGAAATAGGAACACACTGAGTCATCATGAGGAGACTTTG
CAGCATGAGAAGATTCAACTTTAGAGCACAATTTGAATACAGTATATGTCAGGAAACCTCCTTGAA
AAGGCAGTATTCAATACACAGAAGAGAGAGAACCGAGAAGAGAATAACTGTGATTATAATGAATTTGGG
AGAACTTTGTGTGATAGTTCATCCCTCTTGTCCATCAGATATCTCCGTCAAGGGACAATCACTATGAA
TTTAGTGATTGTGAGAAGTCTTATGTGTGAAGTCCACCCCTTTCTAAACCTCATGGGGTATCTATGAAA
CACTATGATTGTGGTGAAGTGGGAATAATTTCAAGGAGGAAATTTGTGTCTGCACACCTTCAGAAAGGT
GATAAAGGAGAGAAAACACTTTGAATGTAATGAATGTGGGAAAGCTTTCTGGGAGAAGTCACATCTCACT
CGACATCAGAGGGTGCACACAGGACAGAAAACCTTTCAATGTAATGAATGTGAAAAGCTTTCTGGGAT
AAGTCAAACCTCACTAAACATCAAAGATCACACACAGGGGAGAAAACCTTTTGAATGCAATGAATGTGGG
AAAGCTTTAGCCATAAGTCAGCCCTCACATTACACCAGAGAACACATACAGGGGAGAAAACCTTCAA
TGTAATGCGTGTGGGAAAACCTTTTCCAGAAATCTGACCTCACTAAACATCAGAGAACACACACAGGG
CTGAAAACCTATGAATGTTATGAATGTGGAAAATCCTTCCGTGTGACTTCGCACCTTAAAGTACACCAG
AGAACTCACACAGGTGAGAAAACCTTTTGAATGTCTTGAGTGTGGGAAATCCTTTAGTGAAGTCAAAAT
CTTACACAGCATCAGAGAATTCACATAGGAGATAAATCTTATGAATGTAATGCATGTGGGAAAACCTTC
TACCACAAGTCATTACTCACCAGGCATCAGATAATTCATACAGGGTGGAAAACCTTATGAATGTTATGAA
TGTGGGAAAACCTTCTGCTTGAAGTCAGACCTCACAGTACATCAGAGAACACACACAGGGCAGAAAACCC
TTTGCATGTCCCGAATGTGGGAAATCTTTAGCCATAAGTCAACCCTCTCTCAACATTATAGAACACAC
ACAGGGGAGAAAACCTACGAATGTCATGAATGTGGAAAATCTTTACAATAAATCATACCTAACTAAA
CATAATAGAACACATACAGGGGAGAAAACCTATGAATGTAATGAATGTGGAAAAGCTTCTACCAGAAG
TCAAACTCACTCAGCATCAGAGAATTCACATAGGGGAGAAAACCTATAAATGTAATGAGTGTGGAAAA
GCTTTCTGCCATAAGTCAGCTCTAATTGTACATCAGAGAACCATACACAAGAAAAGCCCTATAAATGT
AATGAATGTGGAAAATCTTTCTGTGTAATAACAGGACTTATTTCCATGAGAGAAAAGCACACGGGGGAG
AAACCCTATGAATGCAATGAATGTGGGAAATCTTTCAGGCACAAATCATCACTCACAGTACATCACAGG
GCTCACACAGGAGAGAAAATCTTGTCAATGTAATGAATGTGGAAAATCTTTACCGTAAATCGGAACTT
GCTCAACATCAGAGATCACATACAGGGGAAAAGCCCTATGAATGTAACACATGCAGGAAAACCTTTCTCT
CAAAAGTCAAATCTCATTGTACATCAGAGAAGACATATAGGAGAAAACCTTATGAATGAAATGGATATT
AGAAATTTCCAGCCACAAGTCAGCCTCCATAATGCCTCAGAGTATTCACACTGTGGAGAAAAGCCCTGAT
GACATCCTGAATGTTCAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** Sgfl-Mlul

## Plasmid Map:



ACCN: NM\_006954

Insert Size: 2436 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_006954.1](#)

RefSeq Size: 6130 bp

RefSeq ORF: 2436 bp

Locus ID: 7581

UniProt ID: [Q06730](#)

Cytogenetics: 10p11.1

**Protein Families:** Transcription Factors, Transmembrane

**MW:** 94.5 kDa

**Gene Summary:** May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]  
Transcript Variant: This variant (1) represents the predominant transcript and encodes isoform a.