

Product datasheet for **SC311468**

WDR38 (NM_001045476) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	WDR38
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PSI00001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC311468 representing NM_001045476. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAAATTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAACAGCGGGGTCCCAGCCACGCTGGCCGTGCGGAGAGTGAATTCTTCGGCCAGCACGGCGGGGAG
GTCAACTCTTCGCTTCTCCCTGATGGCCAGATGCTGCTCACAGGCTCAGAAGATGGCTGCGTGTAT
GGCTGGGAGACCCGAGTGGGAGCTGCTGTGGAGGCTGGGTGGCCACACAGGCCCGTGAAGTTCTGC
CGCTTCTCCCTGATGGCCACCTCTTCGCCAGCGCCTCCTGTGACTGCACTGTCCGCCTGTGGGATGTG
GCAAGAGCGAAGTGTCTGCGGGTCTGAAGGGTCACCAACGGAGTGTGGAGACGGTCAGCTTCAGCCCT
GACTCGAGACAGCTGGCATCAGGTGGCTGGGACAAGCGGGTGATGCTCTGGGATGTGCAGTCCGGCCAG
ATGCTGCGCCTCTTAGTTGGGCACCGTACTCCATCCAGAGCAGCACTTCTCACCACGGTGAAGTGC
CTGGCCACCGGCTCCTGGGACTCCACCGTACACATCTGGGACCTGCGGATGGTGACCCAGCAGTCTCC
CACCAGGCGTAGAGGGACACAGTGCCAACATCAGCTGCCTGTGCTATTGAGCATCCGGCCTCCTGGCA
TCCGGCTCCTGGGACAAGACCATCCACATCTGGAAGCCACAACCAGCAGCCTGCTTATCCAAGTGAAG
GGCCATGTACCTGGGTGAAGAGCATAGCCTTCTCTCCGACGAGCTGTGGCTGGCCAGCGCCGGCTAT
TCCCGCATGGTCAAAGTCTGGGACTGCAACACAGGAAAGTGCCTTGAGACCCTGAAGGGAGTCTGGAT
GTGGCCACACCTGTGCCTTCAACCCAGATGGGAAATCTTAGTGTCTGGAGCTGCCGATCAGACAGTA
CGTCAAATATCCCGCACGTCAAATCACCAGGGACCTCAAACCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGCCCGGC
```

Restriction Sites:	SgfI-MluI
ACCN:	NM_001045476
Insert Size:	945 bp



OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001045476.2</u>
RefSeq Size:	1325 bp
RefSeq ORF:	945 bp
Locus ID:	401551
UniProt ID:	<u>Q5JTN6</u>
Cytogenetics:	9q33.3
MW:	34.3 kDa