

Product datasheet for SC311469

MPPED1 (NM_001044370) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MPPED1 (NM_001044370) Human Untagged Clone
Tag:	Tag Free
Symbol:	MPPED1
Synonyms:	239AB; C22orf1; FAM1A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC311469 representing NM_001044370. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
 ATGTGGCGCTCTAGTGGGATGCCAGCGTCCTGAAGGCGGAGGCCCTGGCCCTCCTCCCCTGCGGCCTG
 GGCATGGCATTCTCCAGTCCCACGTGATGGCCGCTCGGCGGCACCAGCACAGCCGGCTCATCATCGAG
 GTGGACGAGTACAGCTCCAACCCACCCAGGCCTTCACCTTCTACAACATCAACCAGGGCCGCTTCCAG
 CCACCGCATGTGCAGATGGTGGACCCGGTGCCTCACGATGCCCCAAACCTCCAGGCTACACCCGCTTC
 GTCTGCGTCTCTGATACCACTCGAGGACGACCCCATCCAGATGCCGTACGGCGACGTGCTGATCCAC
 GCTGGGACTTCACTGAGCTGGGCTCCCGAGCGAGGTGAAGAAGTTCAACGAGTGGCTGGGCAGCCTG
 CCCTACGAGTACAAGATCGTGATCGCAGGCAACCACGAGCTGACCTTTGACCAGGAGTTCATGGCCGAC
 CTCATCAAGCAGGACTTTTACTACTTCCCATCTGTGTGCGAAGCTGAAGCCGGAGAACTATGAGAATGTG
 CAGTCGCTGCTGACCAACTGCATCTACCTTCAGGACTCGGAGGTACCGTGCGGGGCTTCCGGATCTAT
 GGCTCCCCATGGCAGCCCTGGTTCTACGGCTGGGCTTCAACCTCCCGCAGGCCAAGCCCTGCTGGAG
 AAATGGAACCTCATTCCGAAGGCGTAGACATCCTGATAACCCATGGACCACCACTGGGCTTCTGGAC
 TGGGTCCCAAGAAGATGCAGCGGTGGGCTGTGTGGAGCTGCTCAACACGGTGCAGAGGCGCGTCCAG
 CCGCGGTTACATGTCTTTGGCCACATCCACGAAGGGTATGGTGTGATGGCAGATGGGACGACCACCTAT
 GTGAATGCGTCCGTATGCACTGTGAACCTACCAGCCGTAACCCGCCCATAGTCATCGACCTCCCACA
 CCGCGAACTCCTGA
 AC**CGGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001044370
Insert Size:	981 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001044370.1</u>
RefSeq Size:	2531 bp
RefSeq ORF:	981 bp
Locus ID:	758
UniProt ID:	<u>O15442</u>
Cytogenetics:	22q13.2
Protein Families:	Transcription Factors
MW:	37.2 kDa
Gene Summary:	May have metallophosphoesterase activity (in vitro).[UniProtKB/Swiss-Prot Function]