

Product datasheet for SC311475

CSNK1G3 (NM_001044722) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CSNK1G3 (NM_001044722) Human Untagged Clone
Tag:	Tag Free
Symbol:	CSNK1G3
Synonyms:	casein kinase 1, gamma 3; OTTHUMP00000159174; OTTHUMP00000159175
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC311475 representing NM_001044722. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAAAATAAAAAGAAAGACAAGGACAAATCAGATGATAGAATGGCAGCAGCTAGTGGTCGATCGGGA
CACAACTCGAGGAACTGGGTCTTCATCGTCTGGAGTTTTAATGGTTGGACCTAACTTTAGAGTTGGA
AAAAAATTTGGATGTGGCAATTTTGAGAAATTACGATTAGGGAAAAATTTATACAAAATGAATATGTG
GCAATTAAGTTGGAGCCCATGAAATCAAGAGCACCACAGCTACATTTGGAATACAGATTCTATAAGCAG
TTAGGATCTGGAGATGGTATACCTCAAGTTTACTATTTTCGGCCCTTGTGGTAAATACAATGCTATGGTG
CTGGAAGTCTGGGACCTAGTTTGAAGACTTGTGGTACTTGTGTGACAGAACATTTTCTCTAAAACA
GTTCTCATGATAGCTATACAAGTATTTCTCGCATGGAATATGTCCATTCAAAGAACTTGATATACAGA
GATGTAACCTGAGAACTTCTTAATAGGACGACCAGGAAACAAAACCCAGCAAGTTATTCACATTATA
GATTTTGGTTTGGCAAAGGAATATATTGATCCGGAGACAAAGAAACACATACCATACAGAGAACAAG
AGCCTTACAGGAACAGCTAGATATATGAGCATAAACACACATTTAGGAAAAGAACAAAGTAGAAGAGAC
GATTTAGAAGCTTTAGGTCATATGTTTCTGAGAGGCAGTCTTCTTGGCAAGGCTTAAAG
GCTGACACATTAAGGAGAGGTATCAGAAAATGGAGATACAAAACGGGCTACACCAATAGAAGTGTTA
TGTGAAAATTTCCAGAAATGGCAACATATCTTCGTTATGTAAGAAGGCTAGATTTTTTTGAAAACCA
GACTATGACTACTTAAGAAAGCTTTTTACTGACTTGTGGTGCAGAAAGGATATATGTTTGTATGAA
TATGACTGGATTGGTAAACAGTTGCCTACTCCAGTGGGTGCAGTTCAGCAAGATCCTGCTGTGCATCA
AACAGAGAAGCACATCAACACAGAGATAAGATGCAACAATCCAAAAACCGGTTGTAAGTTCTACAAT
GGAGAGTTAAACACAGATGACCCACCAGGACGTTCAAATGCACCCATCACAGCCCCTACTGAAGTA
GAAGTGATGGATGAAACCAAGTGTCTGTTTTTTCAAACGAAGGAAAAGGAAAACCATACAGCGCCAC
AAATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001044722

Insert Size: 1248 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001044722.1](#)

RefSeq Size: 4342 bp

RefSeq ORF: 1248 bp

Locus ID: 1456

Cytogenetics: 5q23.2

Protein Families: Druggable Genome, Protein Kinase

Protein Pathways: Hedgehog signaling pathway

MW: 47.9 kDa

Gene Summary: This gene encodes a member of a family of serine/threonine protein kinases that phosphorylate caseins and other acidic proteins. A related protein in the African clawed frog participates in the transmission of Wnt/beta-catenin signaling. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jul 2012]