

Product datasheet for SC311476

ZNF557 (NM_001044388) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF557 (NM_001044388) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF557
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC311476 representing NM_001044388. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
 ATGGCGGCTGTCGTCTGCCCAACTGCCGCTCTCAGCGAGAAGGACACACAGAGGGCGGAGAGCTG
 GTTAATGAGCTCTGAAAAGCTGGCTAAAGGGCTTGGTGACCTTTGAGGATGTGGCCGTGGAGTTCAAC
 CAGGAGGAGTGGGCATTGCTGGACCTGCCAAAGGACACTGTACAGGGACGTGATGCTGGAGAAGTGC
 AGGAACCTGGCCTCACTGGGGAACCAAGTTGATAAACCTAGGCTGATCTCCAGCTGGAGCAAGAAGAT
 AAAGTGATGACAGAAGAGAGAGGAATTCCTCAGGTACCTGTCCAGATGTGGAGAATCCATTTAAAGCC
 AAAGGGTTAACTCCTAAGCTGCATGTTTTTCGAAAAGAACAATCTAGAAATATGAAAATGGAGAGGAAT
 CATCTTGGAGCAACACTCAACGAATGTAATCAGTGTTTTAAAGTCTTCAGCACAAAATCTCCCTTACA
 CGGCACAGGAAGATTCATACTGGAGAAAGACCCTATGGCTGCAGTGAATGTGGGAAATCCTACAGCAGT
 AGATCTTACCTTGCTGTTTCATAAGAGAATCCACAATGGGAGAAAACCTATGAATGCAATGACTGTGGG
 AAAACCTTCAGCAGCAGATCTTACCTTACTGTTTCATAAGAGAATCCACAATGGGAGAAAACCTACGAA
 TGCAGTGACTGTGGGAAAACCTTCAGCAATTCCTCATACCTCAGACCGCACTTGAGAATTCACACTGGA
 GAAAAACCGTACAAATGTAACCAAGTGTTCGTGAGTTCGCACTCAGTCAATCTTCACAAGGCACAAG
 AGAGTTTCATACGGGGAGGGTCATTATGTATGTAATCAGTGTGGAAGGCTTTCGGCAGCAGGTCATCT
 CTTTCTTCGCACTATAGCATTATACAGGGGAGTACCCTTACGAATGCCACGATTGTGGGAGAACCTTC
 AGGAGGAGGTGCAATCTGACACAGCACATAAGAACTCATCTGGAGAAAAACCTACACATGTAATGAG
 TGTGGGAAATCCTTACCAATAGCTTTCTCTTACAATTCACAGGAGAATACATAATGGAGAGAAATCC
 TATGAGTGCAGTGATTGTGGAAATCCTTAAATGTTCTCTCATCCGTTAAGAAACACATGAGAACTCAC
 ACTGGAAAAAACCTATGAATGTAATTATTGCGGGAAATCCTTCACAAGTAACCTACCTTCTGTG
 CATACGAGAATGCATAATAGGCAATGTGA
 ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites: SgfI-MluI


[View online »](#)

ACCN:	NM_001044388
Insert Size:	1272 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001044388.1</u>
RefSeq Size:	5981 bp
RefSeq ORF:	1272 bp
Locus ID:	79230
UniProt ID:	<u>Q8N988</u>
Cytogenetics:	19p13.2
Protein Families:	Transcription Factors
MW:	48.6 kDa
Gene Summary:	<p>May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (3) uses an alternate splice site in the 5' coding region compared to variant 1. The resulting protein (isoform b) is shorter but has the same N- and C-termini compared to isoform a.</p>