

Product datasheet for **SC311480**

SMAP1 (NM_001044305) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SMAP1 (NM_001044305) Human Untagged Clone
Tag:	Tag Free
Symbol:	SMAP1
Synonyms:	SMAP-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC311480 representing NM_001044305.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCGACGCGCTCCTGTCTGGGAGAAGGCTCAGAAGCTGAACGAGCAGCACCAGCTCATCCTATCCAAG
CTTCTGAGGGAGGAGACAACAAGTACTGCGCCGACTGCGAGGCCAAAGTCTCGATGGGCTTCTGG
AATATTGGTGTGTTTATTTGCATCAGATGTGCTGGAATTCATAGAAATCTTGGGGTTCATATATCCAGG
GTCAAATCAGTCAACCTAGACCAATGGACAGCAGAACAGATACAGTGCATGCAAGATATGGGAAATACT
AAAGCAAGACTACTCTATGAAGCCAATCTTCCAGAGAAGTTCGAAGACCACAGACAGATCAAGCAGTG
GAATTTTTTCATCAGAGATAAATATGAAAAGAAGAAATACTACGATAAAAAATGCCATAGCTATTACAAT
ATTTCTCTCTGATGCTCCTCTCAGCCTTTGGTATCCTCTCCTCTCTGCAAGCTGCTGTTGACAAA
AATAAATTGGAGAAAGAAAAGAAAAAAGGAAGAGAAAAAGAGAAAAAGGAGCCAGAAAAGCCG
GCAAAACCACTTACAGCTGAAAAGCTGCAGAAGAAAGATCAGCAACTGGAGCCTAAAAAAGTACCAGC
CCTAAAAAGCTGCGGAGCCACTGTGGATCTTTTAGGACTTGATGGCCCTGCTGTGGCACCAGTGACC
AACGGGAACACAACGGTGCACCCCTGAACGATGATCTGGACATCTTTGGACCGATGATTTCTAATCCC
TTACCTGCAACTGTCATGCCCCAGCTCAGGGACACCCTCTGCACCAGCAGCTGCAACCCTGTCTACA
GTAACATCTGGGGATCTAGATTTATCTACTGAGCAAATACAAAAACAGAAGAAGTGGAAGAAACAA
CTTTCAAAGACTCCATCTTATCTGTATGGCACAGGAACCATCAACAGCAAAGTACTCCTGGTGTA
TTATGGGACCCACAAATATACCATTTACTCACAGCACCAGCTGCATTTACAGGGCTTCCATCGATG
GGCGTGCCTGTGCTGCAGCTCCTGGCCTTATAGGAAATGTGATGGGACAGAGTCCAAGCATGATGGTG
GGCATGCCCATGCCCAATGGGTTTATGGGAAATGCACAACTGGTGTGATGCCACTTCTCAGAAGCTT
GTTGGCCCCAAGGAGGAATGGTGGGACAAATGGGTGCACCCAGAGTAAGTTTGGCCTGCCGCAAGCT
CAGCAGCCCCAGTGGAGCCTCTCACAGATGAATCAGCAGATGGCTGGCATGAGTATCAGTAGTCAACC
CCTACTGCAGTTTTGGCCAGCCCTCCAGCACACAGCAGGATGGTCTGGAAGCTCATCAGGTCAGACT
CTCAGCACAACTGTGGAAATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN: NM_001044305

Insert Size: 1404 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001044305.2</u>
RefSeq Size:	3349 bp
RefSeq ORF:	1404 bp
Locus ID:	60682
UniProt ID:	<u>Q8IYB5</u>
Cytogenetics:	6q13
Protein Pathways:	Endocytosis
MW:	50.4 kDa
Gene Summary:	<p>The protein encoded by this gene is similar to the mouse stromal membrane-associated protein-1. This similarity suggests that this human gene product is also a type II membrane glycoprotein involved in the erythropoietic stimulatory activity of stromal cells. Alternate splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (1) encodes the longest isoform (A).</p>