

Product datasheet for **SC311487**

C19orf28 (MFSD12) (NM_021731) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C19orf28 (MFSD12) (NM_021731) Human Untagged Clone
Tag:	Tag Free
Symbol:	MFSD12
Synonyms:	C19orf28; PP3501
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC311487 representing NM_021731.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)


```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGCCCGGGACCCACGCGCGGAGCGCGCCGTCGCCCGCGCCGCTGTCCCTGGTGGCGCGGCTG
AGCTACGCCGTGGGCCACTTCCCAACGACCTGTGCGCGTCCATGTGGTTCACCTACCTGCTGCTCTAC
CTGCACCTCGGTGCGCGCCTACAGCTCCCGCGCGCGGGGCTGCTGCTGCTGCTGGCCAGGTGGCCGAC
GGGCTGTGCACACCGCTCGTGGGCTACGAGGCCGACCGCGCCGCGAGCTGCTGCGCCCGCTACGGCCCG
CGCAAGGCCTGGCACCTGGTCGGCACCGTCTGCGTCTGCTGTCTTCCCTTTCATCTTACAGCCCTGC
CTGGGCTGTGGGGCGGCCAGCCCGAGTGGGCTGCCCTCCTCTACTACGGCCCGTTCATCGTGATCTTC
CAGTTTGGCTGGGCTCCACACAGATCTCCACCTCAGCCTCATCCGGAGCTCGTACCAACGACCAT
GAGAAGGTGGAGCTCACGGCACTCAGGTATGCGTTCACCGTGGTGGCCAACATCACCGTCTACGGCGCC
GCCTGGCTCCTGTGCACCTGCAGGGCTCGTCGGGGTGGAGCCACCAAGACATCAGCATCAGCGAC
CAGCTGGGGGCCAGGACGTGCCGTGTTCCGGAACCTGTCCCTGCTGGTGGTGGTGTGCGCGCCGTG
TTCTACTGCTATTCCACCTGGGCACCCGGGAGAGGCGCCGCGCCATGCGGAGGAGCCAGGCGAGCAC
ACCCCTGTGGCCCTGCCACGGCCAGCCCTGCTGCTCTGGAAGCACTGGCTCCGGGAGCCGGCT
TTTACCAGGTGGGCATACTGTACATGACCACCAGGCTCATCGTGAACCTGTCCAGACCTACATGGCC
ATGTACCTCACCTACTCGCTCCACCTGCCCAAGAAGTTTATCGCGACCATTCCCTGGTGTGATGACCTC
AGCGGCTTCTTGTCTCCTTCTCATGAAGCCATCAACAAGTGCATTGGGAGGAACATGACCTACTTC
TCAGGCCCTCCTGGTGATCCTGGCCTTGGCCCTGGGTGGCGTGGCGGAGGACTGGGTGTGGCCGTG
TACGCAGCGGCTGTGCTGCTGGTGTGGTGTGCCACCATCCTCGTCACTCGTGGCCATGACGGCC
GACCTCATCGGTCCCAACACGAACAGCGGAGCGTTCGTGTACGGTCCATGAGCTTCTTGGATAAGGTG
GCCAATGGGCTGGCAGTCATGGCCATCCAGAGCCTGCACCCTTGCCCTCAGAGCTCTGCTGCAGGGCC
TGGTGTAGCTTTTACCCTGGGCGATGGTGGTGTGACGGGCGCGTGGGCGTGGCCGTGCCCTGTGT
CTCTGTAGCTCCTGCTGTGGCCGACCCGCTGCGACGCTCTTTTCTTGCCTGGAGAAGAGGGAGGGGA
GAGGACAAGGGCCCTGGCTACTCCTGGATTCTACAGTCTTGTCCAGCCTCCAAGACCCACAAGTCCC
TTCCTCTGGGAAGCCCCCTGGCTGGAGGTGCACCAGGAAGAAGTGGTCTGGGGCTGGCACTAAGCCA
TGGCCAGGGAAGACTGGGGACCCACTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites:

Sgfl-MluI

Plasmid Map:



ACCN:	NM_021731
Insert Size:	1617 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_021731.2</u>
RefSeq Size:	2284 bp
RefSeq ORF:	1617 bp
Locus ID:	126321
Cytogenetics:	19p13.3
Protein Families:	Transmembrane
MW:	58.7 kDa
Gene Summary:	Transporter that mediates the import of cysteine into melanosomes, thereby regulating skin pigmentation (PubMed:33208952). In melanosomes, cysteine import is required both for normal levels of cystine, the oxidized dimer of cysteine, and provide cysteine for the production of the cysteinyl dopas used in pheomelanin synthesis, thereby regulating skin pigmentation (PubMed:33208952). Also catalyzes import of cysteine into lysosomes in non-pigmented cells (PubMed:33208952).[UniProtKB/Swiss-Prot Function]