

## Product datasheet for **SC312075**

### ZNF596 (AK093442) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF596 (AK093442) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF596
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for AK093442, the custom clone sequence may differ by one or more nucleotides ATGAGATCTCATACTCAAGAGGATCCTTTTCTATGCAATGACTTAGGAGAAGATTTCACT CAACATATAGCATTGACTCAAAATGTGATTACCTACATGAGAACGAAACACTTTGTAAAGC AAAAAGTTGGGAAAATCTTCAGTGACTGGTTATCCTTTAATCAACACAAGGAAATTCAC ACCAAATGTAATCATATGGAAGTCATCTATTTGATTATGCCTTTATCCAAAACCTGCC CTTAGACCACACAGTGTGACTCACACTAGAGAGATAACATTGGAATGTCGTGTGTGGG AAAACCTTTAGCAAAAATTCTAATCTTAGGCGACATGAGATGATTCACACTGGAGAGAAA CCACACGGATGCATCTATGTGGGAAAGCCTTTACTCATTTAATTATTGTGAATACC
Restriction Sites:	Please inquire
ACCN:	AK093442
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>



[View online »](#)

RefSeq: [AK093442.1](#), [BAC04166.1](#)  
RefSeq Size: 1847 bp  
RefSeq ORF: 420 bp  
Locus ID: 169270  
Cytogenetics: 8p23.3  
Domains: zf-C2H2  
Protein Families: Druggable Genome, Transcription Factors  
Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]