

## Product datasheet for **SC312131**

### ZKSCAN1 (AK023770) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ZKSCAN1 (AK023770) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZKSCAN1
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for AK023770, the custom clone sequence may differ by one or more nucleotides
Restriction Sites:	Please inquire
ACCN:	AK023770
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
RefSeq:	<u>AK023770.1</u>
RefSeq Size:	5509 bp
RefSeq ORF:	5509 bp



[View online »](#)

<b>Locus ID:</b>	7586
<b>Cytogenetics:</b>	7q22.1
<b>Protein Families:</b>	Transcription Factors
<b>Gene Summary:</b>	<p>This gene encodes a member of the Kruppel C2H2-type zinc-finger family of proteins. This encoded protein may function as a transcription factor that regulates the expression of GABA type-A receptors in the brain. Transcripts from this gene have been shown to form stable and abundant circular RNAs. Elevated expression of this gene has been observed in gastric cancer and the encoded protein may stimulate migration and invasion of human gastric cancer cells. [provided by RefSeq, Oct 2016]</p>