

Product datasheet for **SC312355**

HERC6 (AK092403) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HERC6 (AK092403) Human Untagged Clone
Tag:	Tag Free
Symbol:	HERC6
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for AK092403, the custom clone sequence may differ by one or more nucleotides ATGACTCTGAATGATATAAAAATAATACAAGTTTCCTGTGGACACTACCACTCCCTGGCA TTATCAAAGATAGCCAAGTGTTTTTCGTGGGAAAGAACAGCCATGGGCAGCTGGGCTTG GGGAAGGAGTTCCCCTCCAAGCCAGCCCGCAGAGGGTGAGGTCCCTGGAGGGGATCCCA CTGGCTCAGGTGGCTGCCGGAGGGGCTCACAGCTTTGCCCTGTCTCTGTGGGACTTCG TTTGGCTGGGGAAGTAACAGTGCCGGGCAGCTGGCCCTCAGTGGGCGTAATGTCCCAGTG CAAAGCAACAAGCCTCTCTCAGTCGGTGCCTGAAGAATCTAGGTGTGGTTTATATCAGC TGTGGTGATGCACACACTGCGGTGCTTACCCAGGACGGGAAAGTGTTCACATTTGGAGAC AATCGCTCTGGACAGCTGGGATACAGCCCCACTCCTGAGAAGAGAGGTCCACAACCTGTG GAAAGAATTGATGGCCTAGTTTCGAGATAGATTGTGGAAGAACTCTAAGCATGGATT ATGACTTTCTTT
Restriction Sites:	Please inquire
ACCN:	AK092403
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [AK092403.1](#), [BAC03879.1](#)

RefSeq Size: 2293 bp

RefSeq ORF: 555 bp

Locus ID: 55008

Cytogenetics: 4q22.1

Domains: RCC1

Protein Families: Druggable Genome

Gene Summary: HERC6 belongs to the HERC family of ubiquitin ligases, all of which contain a HECT domain and at least 1 RCC1 (MIM 179710)-like domain (RLD). The 350-amino acid HECT domain is predicted to catalyze the formation of a thioester with ubiquitin before transferring it to a substrate, and the RLD is predicted to act as a guanine nucleotide exchange factor for small G proteins (Hochrainer et al., 2005 [PubMed 15676274]).[supplied by OMIM, Mar 2008]