

Product datasheet for **SC312491**

SFT2C (SFT2D3) (NM_032740) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SFT2C (SFT2D3) (NM_032740) Human Untagged Clone
Tag:	Tag Free
Symbol:	SFT2C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC312491 representing NM_032740. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
ATGGCGGACCTCCACCGCAGCTGCAGGAGTACCTGGCGCAGGGGAAAGCTGGCGGCCCGCGGCCGCG
GAGCCGCTGCTCGCCGCGGAGAAGGCGGAGAGCCCGGGACCGCCGGCGGAGGAGTGGCTGGGCCG
GCGGGCTTGCGCTGGACGTGGCGCGGAGCCCTGCGGAGTCGGCAGCGGCCGCGCTGACGTGCCTCCCG
AGCGTGACGCGCGGGCAGCGGCTGGCGGCGGGCGGGTGCCTGCTGCTGGCTGCACTGTGTTTCGGC
CTAGCCGCGCTCTACGCACCGGTGTTGCTGCTGCGCGCGCGCAAGTTCGCGCTGCTCTGGTCACTGGGC
TCGGCGCTGGCGTTGGCGGAAGCGCGCTGCTGCGGGCGCGCGCGTGGCGACGCTGCTGCGCTGC
GAAGAAGCGCCGTCGCCGCCGCGCTGCTCTACATGGCAGCGCTGGCGGCCACGCTGTTCCCGCGCTG
GGCCTTCGCAGCACGCTGCTCACGGTGTGGGCGCGGGCGCGCAGGTGGCGCGCTGCTGGCCGCGCTG
GTTGGGCTGCTGCCCTGGGCGGCGGCACCGCGCTGCGCCTCGCACTGGGTGCGCTGGGCGCGCGCGC
GGCCTCGCCAAGGTGCTGCCCGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_032740
Insert Size:	648 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_032740.3
RefSeq Size:	2824 bp
RefSeq ORF:	648 bp
Locus ID:	84826
UniProt ID:	Q587I9
Cytogenetics:	2q14.3
Protein Families:	Transmembrane
MW:	21.8 kDa
Gene Summary:	May be involved in fusion of retrograde transport vesicles derived from an endocytic compartment with the Golgi complex.[UniProtKB/Swiss-Prot Function]