

## Product datasheet for **SC312568**

### Nebulette (NEBL) (AK096540) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Nebulette (NEBL) (AK096540) Human Untagged Clone
Tag:	Tag Free
Symbol:	Nebulette
Synonyms:	LASP2; LNEBL
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	<p>&gt;NCBI ORF sequence for AK096540, the custom clone sequence may differ by one or more nucleotides</p> <pre> ATGGCACTCAACATGAACAACACTCAAAGGCTATGAAAAGAAGCCCTATTGTAATGCACGC TACCCGAAGCAGTCCTTCACCACGGTGGCAGATACACCTGAAAATCTTCGCCTGAAGCAG CAAAGTGAATTGCAGAGTCAGGTCAAGTACAAAAGAGATTTTGAAGAAAGCAAAGGGAGG GGCTTCAGCATCGTCACGGACACTCCTGAGCTACAGAGACTGAAGAGGACTCAGGAGCAA ATCAGTAATGTAAAATACCATGAAGATTTTGA AAAAACAAGGGGAGAGGCTTTACTCCC GTCGTGGACGATCCTGTGACAGAGAGAGTGAGGAAGAACACCCAGGTGGTCAGCGATGCT GCCTATAAAGGGGTCCACCCTCACATCGTGGAGATGGACAGGAGACCTGGAATCATTGTT GCACCTGTTCTTCCCGGAGCCTATCAGCAAAGCCATTCCCAAGGCTATGGCTACATGCAC CAGACCAGTGTGTCATCCATGAGATCAATGCAGCATTACCAAATCTAAGGACCTACCGA GCCATGTACGATTACAGTGCCCAGGATGAAGACGAGGTCTCCTTTAGAGACGGCGACTAC ATCGTCAACGTGCAGCCTATTGACGATGGCTGGATGTACGGCACAGTGCAGAGAACAGGG AGAACAGGAATGCTCCAGCGAATTACATTGAGTTTGTTAAT           </pre>
Restriction Sites:	Please inquire
ACCN:	AK096540
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [AK096540.1](#), [BAC04812.1](#)

**RefSeq Size:** 2672 bp

**RefSeq ORF:** 705 bp

**Locus ID:** 10529

**Cytogenetics:** 10p12.31

**Domains:** NEBU, SH3

**Gene Summary:** This gene encodes a nebulin like protein that is abundantly expressed in cardiac muscle. The encoded protein binds actin and interacts with thin filaments and Z-line associated proteins in striated muscle. This protein may be involved in cardiac myofibril assembly. A shorter isoform of this protein termed LIM nebullette is expressed in non-muscle cells and may function as a component of focal adhesion complexes. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Mar 2010]