

Product datasheet for **SC312642**

APRIL (TNFSF13) (NM_003808) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	APRIL (TNFSF13) (NM_003808) Human Untagged Clone
Tag:	Tag Free
Symbol:	APRIL
Synonyms:	APRIL; CD256; TALL-2; TALL2; TNLG7B; TRDL-1; UNQ383/PRO715; ZTNF2
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL4</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene ORF sequence for NM_003808 edited ATGCCAGCCTCATCTCCTTTCTTGCTAGCCCCAAAGGGCCTCCAGGCAACATGGGGGGC CCAGTCAGAGAGCCGGCACTCTCAGTTGCCCTCTGGTTGAGTTGGGGGGCAGCTCTGGGG GCCGTGGCTTGTGCCATGGCTCTGCTGACCCAACAACAGAGCTGCAGAGCCTCAGGAGA GAGGTGAGCCGGCTGCAGAGGACAGGAGGCCCTCCAGAAATGGGGAAGGGTATCCCTGG CAGAGTCTCCCGGAGCAGAGTTCGATGCCCTGGAAGCCTGGGAGAGTGGGGAGAGATCC CGGAAAAGGAGAGCAGTGCTCACCCAAAAACAGAAGAAGCAGCACTCTGCTCCTGCACCTG GTTCCCATTAACGCCACCTCAAGGATGACTCCGATGTGACAGAGGTGATGTGGCAACCA GCTCTTAGGCGTGGGAGAGGCCTACAGGCCAAGGATATGGTGTCCGAATCCAGGATGCT GGAGTTTATCTGCTGTATAGCCAGTCTGTTTCAAGACGTGACTTTACCATGGGTCAG GTGGTGTCTCGAGAAGGCCAAGGAAGGCAGGAGACTCTATTCCGATGTATAAGAAGTATG CCCTCCCACCCGACCGGCCTACAACAGCTGCTATAGCGCAGGTGCTTCCATTTACAC CAAGGGGATATTCTGAGTGTCAATTCCCGGGCAAGGGCGAAACTTAACCTCTCTCCA CATGGAACCTTCTGGGTTTGTGAAACTGTGA



[View online »](#)

5' Read Nucleotide Sequence:	<p>>OriGene 5' read for NM_003808 unedited</p> <pre> NGGTTAGAATTGTAAACGACTCACTATAGGCGGCCGGAATTCGCACCAGGTTCCCTCCTG GGTGTCACTGGCAGCCCTGTCCTTCTAGAGGACTGGAACCTAATTCTCCTGAGGCTGAG GGAGGGTGGAGGGTCTCAAGGCAACGCTGGCCCCACGACGGAGTGCCAGGAGCACTAACA GTACCCTTAGCTTGCTTTCCTCCTCCCTCTTTTATTTTCAAGTTCCTTTTATTTCTC CTTGCGTAACAACCTTCTCCCTTCTGCACCAGTGCCTTACCCTTACCCGCCGCCAC CTCCTTGCTACCCCACTTTGAAACCACAGCTGTTGGCAGGGTCCCAGCTCATGCCAGC CTCATCTCCTTTCTGCTAGCCCCAAAGGGCCTCCAGGCAACATGGGGGGCCAGTCAG AGAGCCGGCACTCTCAGTTGCCCTCTGGTTGAGTTGGGGGGCAGCTCTGGGGGCCGTGGC TTGTGCCATGGCTCTGCTGACCCAACAACAGAGCTGCAGAGCCTCAGGAGAGAGGTGAG CCGGCTGCAGAGGACAGGAGGCCCTCCAGAATGGGGAAGGGTATCCCTGGCAGAGTCT CCCGGAGCAGAGTTCGATGCCCTGGAAGCCTGGGAGAGTGGGGAGAGATCCCGAAAAG GAGAGCAGTGCTACCCAAAAACAGAAGAAGCAGCACTCTGTCCCTGCACTGGTCCCCA TTACCGCCACCTCCAGGATGACTCCGATGTGACAGAGGTGATGTGGCAACCANCTTTA CGCGTGGCAGAAGCCCTCAGGCCCAAGAATATGTGTTCCGATCCAGGATGCTGGNAGTT ATCTGCTGTATACCACGGCCCTGTTCCAGACGTGACTTTCACCTGGGTCCAGTGGTGGCT CGAGAAGCC </pre>
3' Read Nucleotide Sequence:	<p>>OriGene 3' read for NM_003808 unedited</p> <pre> CGTTGCGGTTGTTACCTTGCTAGGACGCGNGNCGAGCACTGTGGNAGCGGGTCACAG GGNATGCCACCCGGGATCTGTTCCAGGAAAAGCTATGACCGCGGCCCAATCTAGAGTCGA GTTTTTTTTTTTTTTTTTTTTTATCATTATTCTATTTTATTTAACTTTAAGCAGGGCTT GATCAGAAAGAAGAGCTCAGAGGAGTGAGGCAACACCAACACCAGGGCGTACGATACACC CGAGGCTCCTCCGACACCAGTTCAGCCCGGAAGGGCTTGAATGGCCTGCCACGCCGGG GCAATGGAGCGCGTGAAAGCCAAACCTTGGGACCCGAAGAGTAATTTGACTGCGAAGGC ACCGTCATCGTCGGGGCACCCCTGTGGACGGGCGCATAATAGCTCCATCTCTGGAGATA GTGATGGAAGGCTGGGGTCCGATAGATGGAGAGAAGGACGTAGGATAAAGGCTTCTGTGGG GTGAAGAGGCATAGGGCACCTCGGAAGAGGTGTGGCATGGGAGAACCAGGGGTTTTCCCT AACAAAGGGCGGCTTCTGGACTTCGCCAACCGCGGGTGCTTTCTCGAGCGCCGTTA CAGGGGGGGGCTGTTCTGGATGCTTCCAGTGGGCCCGGACCAGGTCTGGGGCNATGCGT CGGGCCGGGTCCTCCGCCCTTATTTTATACTGGCGGGAAGAATTCCGGAGCCTCCTTG GGGAACAAGGAAGCGAGGATTGCTGATATCCAAAGTTCTAGCGGCTGGGGCAATAGAA AAAGGCGAAAGATGGGTGAACCGGGGGGCATACCCTAGGAAAGAGCTCTCGGGACTCC CGGACACTGGTCCGCGCCCTTAAAAATCCAAGGTCTTTGGCTGTGCTTCCAGAGTGGGG ACTGACCCTGGTTCGTCAAG </pre>
Restriction Sites:	NotI-NotI
ACCN:	NM_003808
Insert Size:	1500 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	The open reading frame of this TrueClone was fully sequenced and found to be a perfect match to the protein associated to this reference.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_003808.2 , NP_003799.1
RefSeq Size:	2276 bp
RefSeq ORF:	753 bp
Locus ID:	8741
UniProt ID:	O75888
Cytogenetics:	17p13.1
Domains:	TNF
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane
Protein Pathways:	Cytokine-cytokine receptor interaction
Gene Summary:	<p>The protein encoded by this gene is a member of the tumor necrosis factor (TNF) ligand family. This protein is a ligand for TNFRSF17/BCMA, a member of the TNF receptor family. This protein and its receptor are both found to be important for B cell development. In vitro experiments suggested that this protein may be able to induce apoptosis through its interaction with other TNF receptor family proteins such as TNFRSF6/FAS and TNFRSF14/HVEM. Alternative splicing results in multiple transcript variants. Some transcripts that skip the last exon of the upstream gene (TNFSF12) and continue into the second exon of this gene have been identified; such read-through transcripts are contained in GeneID 407977, TNFSF12-TNFSF13. [provided by RefSeq, Oct 2010]</p> <p>Transcript Variant: This variant (alpha) represents the longest transcript and encodes the longest isoform (alpha).</p>