

Product datasheet for SC312658

FGFR1OP2 (NM_015633) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FGFR1OP2 (NM_015633) Human Untagged Clone
Tag:	Tag Free
Symbol:	FGFR1OP2
Synonyms:	HSPC123-like; WIT3.0
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC312658 representing NM_015633. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGTTGCACAATTGAGAAGGCACTTGCCGACGCTAAAGCTCTTGTTGAAAGATTAAGAGATCATGAC
GATGCAGCAGAATCTCTGATTGAGCAAACACAGCTCTCAACAAGCGAGTAGAAGCCATGAAACAGTAT
CAGGAAGAAATCAAGAACTTAATGAAGTCGCGAGACATCGGCCACGGTCCACGTTAGTTATGGGAATC
CAGCAAGAAAACAGACAAATCAGAGAGTTGCAACAAGAAAACAAAGAATTACGTACATCTCTGGAAGAA
CATCAGTCGGCCTTGGAACCTATAATGAGCAAGTACCGAGAACAATGTTTAGATTGCTAATGGCTAGC
AAAAAAGATGATCCGGGTATAATAATGAAGTTAAAGAGCAGCACTCCAAGATTGACATGGTACATCGT
AACAAAGTCCGAAGGATTCTTCTTGATGCATCTCGACACATCCTTGAAGCACCTCAACATGGACTGGAG
AGAAGGCACTTGGAAGCAATCAGAATGAAGTGAAGCACATGTTGACCAGATAACTGAAATGGCAGCA
GTAATGAGGAAAGCCATTGAAATTGACGAGCAACAGGGTTGCAAGGAACAAGAACGAATATTTCAACTT
GAACAAGAAAACAAAGGCTTGAGAGAGATCCTTCAAATAACTCGAGAATCATTTTTGAACCTAAGGAAA
GATGATGCGTCGGAAAGTACTTCTTTGTCAGCATTAGTGACCAACAGTGATTTGAGTCTGAGGAAGAGC
TGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites:	SgfI-MluI
ACCN:	NM_015633
Insert Size:	762 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_015633.2</u>
RefSeq Size:	3030 bp
RefSeq ORF:	762 bp
Locus ID:	26127
UniProt ID:	<u>Q9NVK5</u>
Cytogenetics:	12p11.23
Protein Families:	Druggable Genome
MW:	29.4 kDa
Gene Summary:	<p>May be involved in wound healing pathway.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) encodes the longest isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>