

## Product datasheet for SC312708

### GLT1D1 (NM\_144669) Human Untagged Clone

#### Product data:

|                           |                                                                                         |
|---------------------------|-----------------------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                                     |
| Product Name:             | GLT1D1 (NM_144669) Human Untagged Clone                                                 |
| Tag:                      | Tag Free                                                                                |
| Symbol:                   | GLT1D1                                                                                  |
| Mammalian Cell Selection: | Neomycin                                                                                |
| Vector:                   | pCMV6-Entry (PS100001)                                                                  |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                                    |
| Fully Sequenced ORF:      | >SC312708 representing NM_144669.<br>Blue=Insert sequence Red=Cloning site Green=Tag(s) |

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG  
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC  
 ATGCGGCTCCTGTTCTGGCGGTGCTGCGGCCACACACCGGCAACGCGGTACGGCCAGCGCGTTCCG  
 GCCCATCTAGAGGCTGCAGGGCACGTGTGCGTTTTGAAGGATGCCTTTGACTTTGAAAGCCGATCTGAG  
 ATTGCAAACCTCATCTTGGCTGAGAACTGCGAGGCTGCCCTGGCTCTTCATCTCTATAGGGGAGGCAGG  
 CTTTTCGAAGGCCACCGAATCCCTTTTGGAGTCATCTTGGTGGAAGTGAATGAAGATGCCAAC  
 CAGGCGGAAAAAACACAGTCATGGGCAGAGTTCTTGAGGAAGCCAGATTGCTGTCGCTTTCACAGAG  
 TCAATGAAGGAAATGGCAAGCGCAGTGGTCCAGTGTTCAGGAAGGAGTGAAGCCAAAGTG  
 AAAAGAGCCGCTGGGTACGATTGATTGGAGAGATGCCTCAAGAAGATCTGCACGCGGTGGTGAAGAAT  
 TGCTTCGCGGTGGTGAATAGCTCTGTCTCTGAAGCATGTCAGCTGCAATTTTGGAGGCAATGGATTTA  
 GAAGTACCGGTATTGGCCAGGAACATCCCCGGGAATGCTGCCGTGGTGAAGCATGAAGTCACAGGGCTA  
 CTGTTTTCCAATCCTCAGGAGTTTGTTCATCTGGCAAAGAGACTGGTTAGTGATCCTGCATTAGAAAAG  
 GAAATCGTAGTGAACGGAAGGGAATACGTGAGAATGTATCATTTCATGGCAGGTGGAAGAGACACCTAC  
 CAACAGCTCATCAGGAAGCTGGAAGGAAGCACTGAAGATTGA  
 ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT  
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC

|                    |           |
|--------------------|-----------|
| Restriction Sites: | SgfI-MluI |
| ACCN:              | NM_144669 |
| Insert Size:       | 801 bp    |


[View online »](#)

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).                                                                                                                                                                          |
| <b>OTI Annotation:</b>        | This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.                                                                                                                                                                                                                                                      |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>RefSeq:</b>                | <u><a href="#">NM_144669.2</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>RefSeq Size:</b>           | 2823 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq ORF:</b>            | 801 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Locus ID:</b>              | 144423                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>UniProt ID:</b>            | <u><a href="#">Q96MS3</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Cytogenetics:</b>          | 12q24.33                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Domains:</b>               | Glycos_transf_1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>MW:</b>                    | 29.5 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |