

## Product datasheet for **SC312718**

### STK22C (TSSK3) (BC035354) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	STK22C (TSSK3) (BC035354) Human Untagged Clone
Tag:	Tag Free
Symbol:	TSSK3
Synonyms:	OTTHUMP0000008862; serine/threonine kinase 22C (spermiogenesis associated); spermiogenesis associated 3; SPOGA3; STK22C; STK22D; testis-specific serine/threonine kinase 22C; testis-specific serine kinase 3
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for BC035354, the custom clone sequence may differ by one or more nucleotides ATGGAGGACTTTCTGCTCTCCAATGGGTACCAGCTGGGCAAGACCATTGGGGAAGGGACC TACTCAAAAGTCAAAGAAGCATTTCACAAAAACACCAAGAAAAGTGGCAATTAAGTT ATAGACAAGATGGGAGGGCCAGAAGAGTTTATCCAGAGATTCTCCCTCGGGAGCTCCAA ATCGTCCGTACCCTGGACCACAAGAACATCATCCAGGTGTATGAGATGCTGGAGTCTGCC GACGGGAAAATCTGCCTGGTGATGGAGCTCGCTGAGGGAGGGGATGTCTTTGACTGCGTG CTGAATGGGGGGCCACTGCCTGAAAGCCGGGCCAAGGCCCTTTCCGTCAGATGGTTGAG GCCATCCGCTACTGCCATGGCTGTGGTGTGGCCACCGGGACCTCAAATGTGAGAACGCC TTGTTGCAGGGCTTCAACCTGAAGCTGACTGACTTTGGCTTTGCCAAGGTGTTGCCAAG TCACACCGGGAGCTGAGCCAGACCTTCTGCGGCAGTACAGCCTATGCTGCCCCCGAGGTG CTGCAGGGCATTCCCACGATAGCAAAAAGGTGATGTCTGGAGCATGGGTGTGGTCTCTG TATGTCATGCTCTGTGCCAGCCTACCTTTTGACGACACAGACATCCCAAGATGCTGTGG CAGCAGCAGAAGGGGTGTCTTCCCCTCATCTGAGCATCTCGGCCGATTGCCAGGAC CTGCTCAAGAGGCTCCTGGAACCCGATATGATCCTCCGGCCTCAATTGAAGAAGTTAGT TGGCATCCATGGCTAGCAAGCACT
Restriction Sites:	Please inquire
ACCN:	BC035354
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">BC035354.1</a></u> , <u><a href="#">AAH35354.1</a></u>
<b>RefSeq Size:</b>	1444 bp
<b>Locus ID:</b>	81629
<b>Cytogenetics:</b>	1p35.1
<b>Protein Families:</b>	Druggable Genome, Protein Kinase
<b>Gene Summary:</b>	This gene encodes a kinase expressed exclusively in the testis that is thought to play a role in either germ cell differentiation or mature sperm function. [provided by RefSeq, Jul 2008]